

# SonoMe VET

# **User Manual**

5C, 10L, 14L, 5CB, 10LB H5C10L, H5C, H10L

**Rev 1.02** 2022-06-28

**Veterinary Use Only** 





# **Revision History**

Revision	Date	Reason for Change
1.00	2021-07-05	Initial Release
1.01	2022-03-03	Change Contact Number
1.02	2022-06-28	· [Ch.2.1.1] Wireless Charger excepted
		· [Ch.4.2.1] Add specification of
		recommended wireless charger
		· [Appendix A] modify specifications
		which is related to FCC
		· [Ch.5] Deletion of Ch.5 (Biopsy (Only
		for 14L)

< 1 / 96 > Doc # : BN\_SMV\_EOP\_01



# **Table of Contents**

Revision History	1
Table of Contents	2
Chapter 1. Introduction	5
1.1 Safety Symbols ······	6
1.2 Intended use ·····	7
1.3 Guidelines for Safe Usage ······	7
1.4 Contraindications	
1.5 Precautions & Warnings for Use	7
Chapter 2. Preparing for Use	
2.1 Unpacking	10
2.1.1 H5C10L	11
2.1.2 H5C & H10L	12
2.1.3 10L, 14L, 10LB, 5C & 5CB	13
2.1.4 Application for each probe	14
2.2 How to Install App	14
2.3 Using Probe	15
2.3.1 Visual Inspection	15
2.3.2 Cleaning Probe	15
2.3.3 Booting	15
2.4 Wireless Connection	16
Chapter 3. Using SonoMe VET······	17
3.1 Basic Information & B Mode	17
3.2 Basic B/Mode ······	19
3.2.1 B/M Mode	19
3.2.2 Color Doppler Imaging Mode	20
3.2.3 Power Doppler Imaging Mode	21
3.2.4 PW(Pulse Wave) Spectral Doppler Mode	22

Doc #: BN\_SMV\_EOP\_01



	3.3 Animal patient Information	23
	3.4 Measurement ······	24
	3.4.1 Basic Measurement	24
	3.4.2 Heart Rate Measurement ·····	27
	3.5 Report	28
	3.6 Storage & Review	29
	3.6.1 Image Storage ······	29
	3.6.2 Video Storage ······	29
	3.6.3 Review image & Video	29
	3.7 Change Wi-Fi Channel ·····	30
Cł	napter 4. Maintenance ······	31
	4.1 Charging Probe	31
	4.1.1 USB Charging	31
	4.1.2 Wireless Charging	31
	4.2 Cleaning and Disinfecting Probes	32
	4.3 Storage	33
	4.4 Troubleshooting	33
	4.5 Disposal ······	34
	4.6 Maintenance & Protection	34
Cł	napter 5. Safety ······	36
	5.1 Important Safety Considerations	366
	5.1.1 Electrical Safety	
	5.1.2 Mechanical Safety	
	5.1.3 Probe Safety Precautions	
	5.1.4 Cyber Security ·····	388
	5.2 Acoustic Power Principle	3939
	5.2.1 Biological Safety	39
	5.2.2 Thermal Index (TI) & Mechanical Index (MI)	39
	5.2.3 Acoustic Power Statement	400
	5.2.4 Acoustic Power Setting	422

Doc # : BN\_SMV\_EOP\_01



Appendix B. Acoustic Output Data	59
Appendix A. Specification ······	4949
5.3.3 Recommended Separation Distance	477
5.3.2 Electromagnetic Immunity	455
5.3.1 Electromagnetic Emission	43
5.3 Electromagnetic Compatibilities	433
5.2.5 ALARA	433



# **Chapter 1. Introduction**

SonoMe VET is the new generation instruments for ultrasonography with the outstanding feature of being wireless. Different from traditional ultrasound systems with a cable connecting the probe to the main unit, no cable appears at the end of the probe of the system.

SonoMe VET is highly integrated with ultrasound image processing, power management and a wireless signal provider to be connected by the main units. The main units differ from traditional devices by using a Tablet PC or Mobile Phone which supports iOS System. The probe acts as a Wi-Fi Access Point and can be connected to a Tablet PC or Mobile Phone.

This manual is intended to provide a thorough overview of the SonoMe VET and should be carefully read before starting to operate the device.



[Figure 1-1] SonoMe VET App & Wireless Probes

< 5 / 92 > Doc # : BN\_SMV\_EOP\_01



# 1.1 Safety Symbols

Symbols	Description	Symbols	Description
<u>^</u>	WARNING: Warnings to prevent a serious accident, or damage to property	<b>₩</b>	Manufacturing date
	Follow the User Manual	SN	Serial Number
⅓	Type BF applied part (Classification based on degree of protection against electric hazard)	H	Direct Current(DC) Voltage source
<b>▲</b> IPX1	Dripping-proof device: Protected against vertically falling water drops	À	CAUTION: Precautions to prevent a minor accident or damage to property
(((•)))	Non-ionizing electromagnetic radiation	***	The manufacturer's name and address are provided.
Z	The waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately	<b>‡</b>	Recycling: Dispose of properly in accordance with all state, province, and country regulations
<u> </u>	This side up: For the duration of shipping/delivery, the carton should face upright	图	Use no hooks: Absolutely no hand hooks should be attached to pull the parcel
	Fragile: Handle with care		Keep away from rain
EC REP	The Authorized European Representative's address is provided	F©	Federal Communications Commission logo

< 6 / 92 > Doc # : BN\_SMV\_EOP\_01



#### 1.2 Intended use

SonoMe VET is designed to obtain ultrasound images, measure and analyze anatomical structures and fluids from animals like cats, dogs and horses

### 1.3 Guidelines for Safe Usage

SonoMe VET is indicated for examining animals. This product is intended for use by, or by the order of, and under the supervision of, a licensed veterinarian who is qualified for direct use of medical devices.

### 1.4 Contraindications

- This product must not be used for ophthalmological applications, or any other use that involves the ultrasound beam passing through the eyeball.
- Do not use this product in the following situations. Doing so may produce images with inaccurate results:
  - Animal patients who have had surgery, which may have changed the composition
    of the examining tissue (for example, a mastectomy), as this could skew or alter
    the measured density.
  - Animal patients whose bodies contain foreign artifacts (for example, implants), in the examining tissue
  - Endo-cavitary use: Defined as introducing a scanner within a (body) cavity or organ.
     For example, atrium, esophagus, rectum, or vagina.

# 1.5 Precautions & Warnings for Use



#### **Precautions:**

- Keep this manual near the product and refer to it when using the product. Bionet is not responsible for any product issues caused by the user's careless operation or negligence in maintenance.
- The product must work in a clean environment, should avoid areas where moisture,



direct sunlight, or a heater is near the product; that are frequently exposed to vibration; that are too dusty or are not properly ventilated; or where chemicals or gases are present. Do not place anything on top of the instrument.

- The product shall be operated in undisturbed conditions to avoid data transmission interruption.
- When there is wireless channel congestion, switch the channel (Refer to Section 3.6 SETTINGS), and then restart the probe.
- Prescription Use. This product must only be used by persons who have sufficient knowledge of and/or a qualification in clinical pathology. Unqualified operators are prohibited from using the product.
- The product shall be repaired only by professionals recognized by Bionet.
- The product does not have shelf life. Its expected use life is 10 years. After 10 years, though the product still works normally, it is recommended to have it checked by Bionet.
- Useless components must be disposed of in compliance with local regulations.
- Be careful when holding the device, for the device is handheld, it may fall.
- The words "Insufficient Storage Space" will appear on the interface to remind the user to clean up space when storage space of the mobile device will be insufficient.



#### Warnings:

- Never use the product in the presence of flammable or anesthetic gas. Doing so may cause an explosion.
- Do not allow the interior of the product to be exposed to, or immersed in, liquid. Otherwise, fire, electric shock, injury, or damage to the product may occur.
- Changes or Modifications not expressly approved by the party responsible could void the user's authority to operate this device.
- If the product does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
  - Reorient or relocate the receiving antenna.
  - Increase the separation between the product and receiver.
  - Connect the product into an outlet on a circuit different from that to which the receiver is connected.



- Consult the dealer or an experienced radio/TV technician for help.
- The device has been evaluated to meet general RF exposure requirement.
- You may lose user settings or Animal patient information files because of physical shocks to the product or internal errors. Therefore, you should backup this information on a regular basis.

< 9 / 92 > Doc # : BN\_SMV\_EOP\_01



# **Chapter 2. Preparing for Use**

FOR YOUR PROTECTION, please read safety instructions completely before turning the power on.

#### **CAUTIONS**

- Too high ultrasonic intensity and / or long exposure time may cause injury.
- Please do not apply the probe of this product to a scope not covered in this manual.

# 2.1 Unpacking

The product is carefully packed to prevent damage during shipment. Before unpacking, please note any visible damage to the outside of the shipping containers.

Items should be checked whether all ordered items have been received. The following table lists the items which should be received with each product.

[Table 2-1] Item List

List	H5C10L	H5C, H10L	10L,14L, 10LB, 5C, 5CB	Remarks
Probe	0	0	0	
USB Cable (1m)		0	0	Micro 5pin
Wireless charger + USB cable	0			USB-C
Quick Guide	0	0	0	
Pouch	0	0	0	
Strap	-	-	0	

Check the items for damage or defects (For example, cracks, broken parts, liquid leaks, sharp edges, etc.). If there is any damage or defects, stop using the probe immediately and contact Bionet's Customer Service Department.

> Doc #: BN\_SMV\_EOP\_01 < 10 / 92 >



## 2.1.1 H5C10L



[Figure 2-1] H5C10L

#	Items	Descriptions
1	Battery Level Indicator	Indicate Battery Level
2	Wireless Signal Status	Show Wireless Signal Status
3	Convex On/Off Indicator	Indicate Convex Array's Activation Status
4	Power On/Off & Freeze	<ul><li>Power ON/OFF: Press button for about 5 sec.</li><li>Freeze: Press button for about 1 sec.</li></ul>
5	Linear On/Off Indicator	Indicate Linear Array's Activation Status

- Switch Probe type: Press power button for about 3 seconds
  - If Convex Array is On, #3 indicator is On.
  - If Linear Array is On, #5 indicator is On.
- Only H5C10L is available for wireless charging



### 2.1.2 H5C & H10L



[Figure 2-2] H5C & H10L

#	Items	Descriptions
1	Power On/Off & Freeze	· Power ON/OFF : Press button for about 5 sec.
		· Freeze : Press button for about 1 sec.
2	Depth Button	Control Depth
3	Gain – Button	Reduce Gain
4	Gain + Button	Increase Gain
5	Wireless Signal Status	Show Wireless Signal Status
6	Battery Level Indicator	Indicate Battery Level



# 2.1.3 10L, 14L, 10LB, 5C & 5CB



[Figure 2-3] 10L, 14L, 10LB, 5C & 5CB

#	Items	Descriptions
1	Power On/Off & Freeze	· Power ON/OFF : Press button for about 5 sec.
I	Fower On/On & Freeze	· Freeze : Press button for about 1 sec.
2	Wireless Signal Status	Show Wireless Signal Status
3	Battery Level Indicator	Indicate Battery Level



# 2.1.4 Application for each probe

Probe	Application	
H5C10L	Convov	Abdomen, Gynecology, Obstetrics, Cardiac, Urology,
	Convex	Kidney
	Linaan	Thyroid, Small Parts, Pediatrics, Vascular, Carotid, Breast,
	Linear	MSK, Nerve
H5C, 5C, 5CB	Abdomen, Gynecology, Obstetrics, Cardiac, Urology, Kidney	
H10L, 10L,	Thyroid, Small Parts, Pediatrics, Vascular, Carotid, Breast, MSK, Nerve	
10LB, 14L		

# 2.2 How to Install App

Download the "SonoMe VET" app via Apple App Store depending on user Tablet PC or Mobile Phone OS.

< 14 / 92 > Doc # : BN\_SMV\_EOP\_01



### 2.3 Using Probe

### 2.3.1 Visual Inspection

Check the probe for cracks, broken parts, liquid leaks, sharp edges, etc.

#### Warning

Using damaged probes may result in electric shocks and other hazards to the Animal patients and/or users. If damaged, stop using them immediately and contact Bionet's Customer Service Department.

### 2.3.2 Cleaning Probe

All probes must be cleaned and disinfected after each use. Refer to "Chapter 4.2. Cleaning and disinfection".

#### Warning

Probes that have not been cleaned or disinfected may cause bacterial and viral infections

### 2.3.3 Booting

Please check the following before scanning

1. The probe should not be abnormally heated during use. If the temperature of the probe is much higher than body temperature, or the surface temperature of the probe exceeds 40 °C (104°F), the probe is stopped

#### Warning

Do not use the probe which is heated abnormally on an Animal patient, as it may burn Animal patient's skin

2. After turning on the power, check whether the functions of ultrasound image, app, button, etc. are working properly.



#### Warning

If there are any of the above problems, there may be a problem with the product, so please contact Bionet's Customer Service Department.

- The Indicators of Wireless Signal Status and Battery Level in the probe are not visible till the probe is turned on.
- Press the power button for a second and turn on. After turning on, you can check the status of battery level. There are 4-Steps to the battery level. (Refer to chapter 4 for how to charge)
- The wireless signal indicator flashes after a few seconds after turning on the probe. This lets you know that the probe is ready to connect with your tablet PC or mobile phone.
- The probe can be turned off by holding the power button down for about 5 seconds. When the probe is turned off, all indicators (wireless signal status and battery level) are turned off, too.

### 2.4 Wireless Connection

If the probe is ready for wireless connection as described previously, turn on Wi-Fi (if not turned on) in the settings of your tablet PC or mobile phone, and find the SSID of the probe.

The SSID is something like "SX-OXXXXXXXXOOO (X = English, O=Number)".

Before connecting, check that the last 6-digits of the 15-digit ID written on the back of the probe match the last 6-digits of the SSID.

When connecting the probe for the first time, you need to enter a password, which is an 11-digit ID written on the back of the probe.

After Wi-Fi is connected, launch the SonoMe VET app. When the app confirms the connection with the probe, the wireless signal indicator on the probe stops blinking. When the connection is complete, start to Scan.



# Chapter 3. Using SonoMe VET

### 3.1 Basic Information & B Mode



[Figure 3-1] B Mode

#### Information:

- 1. Gain +: Increase the brightness of the applicable Operation mode
- 2. Gain -: Decrease the brightness of the applicable Operation mode
- 3. Depth: Adjusts the scanning depth of an image. The range varies depending on the probe in use
- 4. Focus: Adjusts the focus of an image
- 5. Dynamic: Adjusts contrast by changing the ratio of the minimum and maximum values of input signals. The higher the value, the smoother the displayed image.
- 6. Harmonic: Optimize image using various frequency
- 7. Reject: Reduces noise or echoes from an image to make the image clearer
- 8. Image mode

< 17 / 92 > Doc # : BN SMV EOP 01



- 8 1 . 5CB & 10LB: B Mode, B/M Mode
- 8 2 . 5C, 10L, 14L, H5C10L, H5C, H10L
  - : B Mode, B/M Mode, Color Mode, PDI Mode, PW Mode
- 9. Animal patient: Displays the Animal patient Information screen, which is used for selecting an Animal patient ID from the list or entering new Animal patient information
- 10. Freeze: Pauses an image being scanned or reactivates a paused image
- 11. Play: Play saved or stopped video
- 12. Measure: Starts to measure distance, circumference, area, and volume.
- 13. Annotation: Enter a text on the image
- 14. Biopsy: Biopsy function is currently not supported.
- 15. Clear: Deletes text and measurement result displayed on an image.
- 16. (Image) Save: Save an image
- 17. (Video) Save: Save a Video
- 18. Setting: Hide patient and scanning information from screen. Selecting Wi-Fi channel. DICOM set-up. Cine Loop saving options.
- 19. TGC: Adjust the Gain based on the depth (When you click that area, you can find a TGC menu)
- 20. Preset (Application): Select the preset for scanning.

< 18 / 92 > Doc # : BN\_SMV\_EOP\_01



### 3.2 Basic Mode

### 3.2.1 B/M Mode



[Figure 3-2] B/M Mode


#### Information:

1. M Line: The M Line indicates the relative position of the B/M Mode image in the 2D image. You can move the M Line to change the observation area. Use the finger to move the M Line to the right or left.

Doc # : BN\_SMV\_EOP\_01



# 3.2.2 Color Doppler Imaging Mode



[Figure 3-3] Color Doppler Imaging Mode


#### Information:

- 1. Steer: Adjust ultrasound beam angle to minimize loss of color information
- 2. Gain +: Increase the color brightness
- 3. Gain -: Decrease the color brightness
- 4. Move, Zoom: Use the finger to change the location and size of ROI(Region of Interest) box

< 20 / 92 > Doc # : BN\_SMV\_EOP\_01



### 3.2.3 Power Doppler Imaging Mode




[Figure 3-4] Power Doppler Imaging Mode

#### Information:

- 1. Steer: Adjust ultrasound beam angle to minimize loss of color information
- 2. Gain +: Increase the color brightness
- 3. Gain -: Decrease the color brightness
- 4. Move, Zoom: Use the finger to change the location and size of ROI(Region of Interest) box



### 3.2.4 PW(Pulse Wave) Spectral Doppler Mode



[Figure 3-5] PW Spectral Doppler Mode


#### Information:

- 1. Steer: Changes the direction specified in Sample Volume
- 2. Gain +: Increase the pulse brightness
- 3. Gain -: Decrease the pulse brightness
- 4. Volume: Adjust the size of Sample Volume
- 5. Angle: Adjust the angle of Sample Volume




### 3.3 Animal patient Information

When you press "ANIMAL PATIENT( )" button, you can see the Animal patient Information Menu as Figure 3-6.



[Figure 3-6] Animal patient Information Menu



[Figure 3-7] Edit Animal patient Information

When you complete the Animal patient' information, Press "OK" button for save.


< 23 / 92 > Doc # : BN\_SMV\_EOP\_01



#### 3.4 Measurement

#### 3.4.1 Basic Measurement

Press "Measure" button at Freeze in B mode, you can see 10 measurement items (When scrolling, the rest of the measurements are also visible.) as Figure 3-8. The available measurement methods vary depending on the current Operation mode and Application.



[Figure 3-8] Measurement

#### **■ LENGTH**

Choose "LENGTH". Use the Finger and Set to specify both end points of the measurement area. Specify both end points and then the distance between them will be measured automatically.

If the color of both end points is green, you can move or change but if white, the measurement is complete and cannot be changed.

< **24** / **92** > Doc # : BN\_SMV\_EOP\_01





[Figure 3-9] Measure the LENGTH

#### AREA/CIRCUMFERENCE

Choose "AREA/CIRCUM". Using the Finger, set to specify 3 points of the measurement area and then you can see elliptical trajectory. When you click one of 3 points, you can adjust the size and location and the area and circumstance will be measured automatically.

If the color of 3 points is green, you can move or change but if white, the measurement is complete and cannot be changed.

#### ■ ANGLE

Choose "ANGLE". Draw a straight line at first(Please refer to "LENGTH"). Specify 3<sup>rd</sup> point and then measure the angle formed by the points automatically. When you click one of 3 points, you can adjust the location and the area and circumstance will be measured automatically.

If the color of 3 points is green, you can move or change but if white, the measurement is complete and cannot be changed.

#### ■ TRACE

Measure the area of an irregularly shaped object. Choose "TRACE". Set to specify the




starting point for tracing on the contour of the area to be measured. Trace the curve so that the measurement point returns to the start point. And then you can see the area will be measured automatically

#### **Note**

Measurement for GA (CRL), GA (BPD), GA (GS), GA (FL), GA (HC), GA (AC) is ONLY available to "obstetrics" application

The above measurement functions can be fine-tuned using the on-screen virtual trackball. The virtual trackball(please refer to bottom-right corner in Figure 3-9) can be fine-tuned according to the direction of the measurement point.

Up to 4 measurements can be compared at the same time. Delete all measurements using "Clear(X)" button. When you want to delete one specific measurement, select that measurement and then you can see green column as Figure 3-10. After pressing " $\otimes$ " button, you can delete.




[Figure 3-10] Delete specific measurement

< 26 / 92 > Doc # : BN\_SMV\_EOP\_01



#### 3.4.2 Heart Rate Measurement




[Figure 3-11] Heart Rate Measurement

At Freeze in B/M Mode or PW mode, Select "Heart Rate" from the Measurement menu. A green bar appears, with which you can specify the range. To set the measurement range, use the Finger to place the bar at the desired location (In B/M mode, there are 3 bars.; In PW mode, there are 5 bars). The system will automatically measure the heart rate within the measurement range.



### 3.5 Report

When you press "ANIMAL PATIENT( )" button, you can see a "Report" button at bottom-right in Animal patient Information Menu as Figure 3-12.



[Figure 3-12] Animal patient Information Menu

If you select "Report" at the bottom-right of the Animal patient information, a report window appears as shown in Figure 3-13 and click the "Obs(Observations)" box to allow the user to input observations or opinions. You can download the report by selecting the download icon

"  $\underline{\hspace{1cm}}$  " at the bottom-right of the page.

The report is automatically saved to your phone or tablet's photo album.



[Figure 3-13] Report



### 3.6 Storage & Review

### 3.6.1 Image Storage

As shown in Figure 3-14, if you select "Save Image ( )" at the bottom right of the App screen, the image currently displayed on the screen is saved. The saved image is automatically saved to your phone or tablet's photo album.



[Figure 3-14] Complete to save image

### 3.6.2 Video Storage

As shown in Figure 3-14, if you select "Save Video ( )" at the bottom-right of the App screen, the saved video is automatically saved in the photo album of your phone or tablet PC within 100 seconds.

### 3.6.3 Review image & Video


Open the photo album on your phone or tablet and view the saved images and videos.

< 29 / 92 > Doc # : BN\_SMV\_EOP\_01



### 3.7 Change Wi-Fi Channel

In a congested Wi-Fi environment, the user can choose a different Wi-Fi channel for the probe. Press the "Setting ( )" button, then select the appropriate channel selection by referring to the signal channel selection list (see Figure 3-15). Restart the probe after 2 seconds and reconnect with the mobile phone or tablet according to chapter 3.1.



[Figure 3-15] Wi-Fi Chanel Setting

< **30 / 92 >** Doc # : BN\_SMV\_EOP\_01



# Chapter 4. Maintenance

### 4.1 Charging Probe

### 4.1.1 USB Charging

When the battery is low, the probe needs to be charged. Remove the rubber cover from the tip of the probe and then plug one end of the USB cable into the USB charging port in the probe as shown in Figure 4-1. The battery level indicator blinks while charging.

The battery is fully charged when both the battery level indicator and the wireless link indicator are not blinking. When charging is complete, remove the USB cable and then put the rubber cover into the USB charging port to avoid liquid and device damage.



[Figure 4-1] Probe Charging via USB Cable

### 4.1.2 Wireless Charging

- Plug one end of the USB cable into the wireless charging pad. Place the probe on the wireless charging pad. The probe lights up and starts charging after 1~2 seconds.
- Ensure the probe is centered on the wireless charging pad. For efficient charging, place the probe, place the probe within 7mm left and right from the center of the wireless charging pad. If it is significantly out of the range, it may be unstable or may not charge.

< **31** / **92** > Doc # : BN SMV EOP 01



Wireless Charger Information

Input Power	5V/2A, 9V/2A
<b>Output Power</b>	5W / 10W



[Figure 4-2] Wireless Charging with wireless charging pad

#### Warning

It should not be used if the adapter power supply voltage is outside the adaptive range of the appliance (Normal Adapter's output voltage :  $5V \pm 0.5V$ ).

Check the condition of power cable at first. If the power cable is damaged or broken, replace with a new one immediately.

### 4.2 Cleaning and Disinfecting Probes

All probes must be cleaned and disinfected before and after each use. Use a soft cloth or appropriate cleaning sheet, lightly dampened with Isopropyl alcohol (or proper cleaner), to remove any foreign substances left on the probe, the edges, corners, and curved parts of the probes. Dry the probe with a clean, soft cloth. Alternatively, dampen a soft cloth in a glutaraldehyde-based hospital disinfectant solution such as Cidex. Wipe the probe with a dampened cloth.

If still wet or left a stain, wipe with a clean water-dampened cloth. Before use, dry the probe completely with clean, soft cloth.

> < 32 / 92 > Doc #: BN\_SMV\_EOP\_01



## 4.3 Storage

When not in use, it is recommended to put the probe in its case. During storage, the probe must be protected from extreme operating temperatures

## 4.4 Troubleshooting

Inspection: Check whether connection between probe and host(smart phone or tablet) is OK or not

#### Error handling:

#	Symptoms	Solutions
1	No response after pressing power	Check the status of charging and
'	button	power supply
	Wi-Fi is NOT connected between the	Check whether Wi-Fi signal channel is
2	probe and host (smart phone or	ready or not.
	tablet)	Check whether Wi-Fi password is
	tablet)	correct or not.
		Check whether other equipment
3	Displayed on the screen with noise	started which cause electromagnetic
3	and interference like snow	interference. If so, shut down the
		device or get far from the device.
4	The brightness of image is not	Adjust the brightness and gain
	enough	Adjust the brightness and gain
		Check whether the USB interface,
5	Charging does NOT work	circuit and electrical outlet are
		damaged or not.

< **33** / **92** > Doc # : BN\_SMV\_EOP\_01



### 4.5 Disposal

#### Warning

Do not dispose of the product carelessly.

- Battery recycling must meet local requirements.
- Waste must be disposed of in compliance with local regulations.

#### 4.6 Maintenance & Protection

- 1. This product usage and storage conditions shall comply with the environmental conditions of Appendix A in this manual.
- 2. The product power supply shall be in accordance with Appendix A in this manual.
- 3. If the product is not used for a long time, charge it at least 2 times a week for at least 1 hour.
- 4. Do not tamper with the probe cover to clean or disassemble the components inside the probe.
- 5. Clean the cover with alcohol-cotton after turning off.
- 6. This product should not be started up and shut down too often. If you need to restart, wait at least 1 minute for booting after shutdown.
- 7. If any problem is happened, stop using the device and contact Bionet's Customer Service department.
- 8. Since the probe is vulnerable to external impact, collision or falling is prohibited.
- 9. If you want to pause during scanning, use the "Freeze" function. In freeze status, the device can be kept on standby for a long time.
- 10. Only use the appropriate ultrasound gels for medical standard.
- 11. The probe is waterproof but, do not immerse the probe in conductive liquid to minimize corruption.
- 12. Check for cracks regularly to avoid liquid immersion and damage of internal components.
- 13. Refer to chapter 4.2 about how to clean and disinfect.
- 14. To maintain the performance and safety of the system, regular electrical and



mechanical safety inspections of the system must be performed by professional technicians in less than six months.

15. Excessive shocks, such as a fall, may seriously damage the product. In this case, contact Bionet's Customer Service Department for maintenance and repair.



## **Chapter 5. Safety**

To prevent damage of the equipment or injury to yourself or others, read the following safety instruction before using SonoMe VET.

## **5.1 Important Safety Considerations**

Read the following safety instruction before using SonoMe VET.



- Do not use for any purpose other than the intended use. Otherwise, the system may be damaged or serious injury may result.
- This device is for diagnostic use only and is not intended for therapeutic use.

Always keep this manual with the system. Periodically review operating procedures and safety precautions.

## **5.1.1 Electrical Safety**

- The biocompatibility has been verified so that it does not harm users or Animal patients under normal circumstances
- Never attempt to modify the product in any way.
- The battery is a consumable and will lose performance over time. If the battery life becomes less than 2 hours, it is time for a replacement. When you want to purchase or replace a battery, please contact Customer Service Department.
- Warning: Class I The product must be connected to a power supply with a protective ground to prevent electric shock.
- Do not pour liquids on the system surface. Liquids seeping into electrical circuits can cause excessive leakage currents or system failure. If water is unavoidably poured into the system, immediately stop using the system and contact a Customer Service Department.
- Use only probes provided by the manufacturer. Otherwise, the system will not work and the worst accidents such as fire may occur.



- Check the probe for cracks, broken parts, liquid leaks, sharp edges, etc. If there is any damage, stop using the probe immediately and contact Customer Service Department.



#### Warnings:

- This product is intended for use by, or by the order of, and under the supervision of, a licensed veterinarian who is qualified for direct use of medical devices.
- This system may only be maintained by persons approved or trained by Bionet.
- Always use application-specific probes to obtain the best quality images.
- Never use the product in the presence of flammable or anesthetic gas. Doing so may cause an explosion.
- Do not use the system simultaneously with other equipment such as electric knives, defibrillators, and other high-frequency therapy equipment. Otherwise, there is a risk of electric shock
- Keep the system dry and do not transport to sites with rapid environment changes to avoid short circuits caused by condensation or water droplets.
- The product must be connected to a power supply with a protective ground to prevent electric shock.
- After reading the manual thoroughly, set and control the acoustic power level.

## 5.1.2 Mechanical Safety



#### **Warnings:**

- Be careful when holding the device, for it is handheld, it may fall.
- Do not use shell cracking equipment.

## **5.1.3 Probe Safety Precautions**

- Please only use the ultrasound gels sold legally. Always follow the instructions in the User Manual. To prevent contamination, use and take care of ultrasound gel properly.



#### Warnings:

- Disconnect the probe from the system after freezing the image or turning the system down. Otherwise, the system or probe may be damaged.

< **37 / 92 >** Doc # : BN SMV EOP 01





### Warnings:

- Check the probe for cracks, broken parts, liquid leaks, sharp edges, etc. If there is any damage, stop using the probe immediately. Using damaged probes may result in electric shocks and other hazards to the Animal patients and/or users.
- After disinfecting the accessory, it is necessary to flush the chemicals out of the accessory. Residual chemicals or gases left behind can not only damage the accessory but can also be harmful to animal health.

## **5.1.4 Cyber Security**

- Backup the database regularly to avoid loss or damage of the database.
- The probe can connect to Tablet PC or Smart Phone via local wireless network. The software itself cannot connect to external network. The network which is connected to the software is a local wireless network connected by the probe.
- If there are any software bugs during use, users can give feedback through Bionet e-mail account(sales@ebionet.com) and Bionet will analyze and fix the bug, and if an update is required, the user will be notified of it by posting on the website or e-mail.

< 38 / 92 > Doc # : BN\_SMV\_EOP\_01



## **5.2 Acoustic Power Principle**



#### Warnings:

- Guidance for the use of diagnostic ultrasound is defined by the "As Low As Reasonably Achievable (ALARA)" principle. By keeping ultrasound exposure as low as possible while obtaining diagnostic images, users can minimize ultrasonic bioeffects.
- The operator should notice the effect of the heat on the Animal patient body when the exam is performed around the bones and the nearby soft tissues which can transform the ultrasound energy to heat energy. Take special care to the fetal patient whose bones are growing.

## 5.2.1 Biological Safety

Ultrasound waves may have damaging effects on cells and, therefore, may be harmful to the Animal patient. If there is no medical benefit, minimize the exposure time and maintain the ultrasound wave output level at a low setting. Please refer to the ALARA principle.

## 5.2.2 Thermal Index (TI) & Mechanical Index (MI)

The system output display comprises two basic indices: a thermal index (TI) and a mechanical index (MI). These shows on the Right-Top side on the app.

#### Meaning of MI & TI

Mechanical bioeffects are threshold phenomena that occur when a certain level of output is exceeded. The threshold level varies, however, with the type of tissue. The MI accounts for these two factors. The higher the MI value, the greater the likelihood of mechanical bioeffects occurring. However, there is no specific MI value that means that a mechanical effect will occur. The MI should be used as a guide for implementing the ALARA principle.

The TI informs the user of the potential for temperature increase occurring at the body surface, within body tissue, or at the point of focus of the ultrasound beam on bone. The TI is an estimate of the temperature increase in specific body tissues. The actual amount of any temperature rise is influenced by factors such as tissue type, vascularity, and mode of operation. The TI should be used as a guide for implementing the ALARA principle. Depending on the examination and type of tissue



involved, TI could be one of three types.

- The soft tissue thermal index (TIs): the user about the potential for heating within soft homogeneous tissue.
- The bone thermal index (TIb): the user about potential heating at or near the focus after the ultrasound beam has passed through soft tissue or fluid, such as the skeletal structure of a 2-3 month-old fetus.
- The cranial bone thermal index (Tlc): the user about the potential heating of bone at or near the surface, for example, the cranial bone.

#### - Mechanical and Thermal Indices Display Precision and Accuracy

TI and MI values are displayed in real time on the screen. The operator should observe these index values during examinations and ensure that exposure time and output values are maintained at minimum amounts needed for effective diagnosis. The Mechanical and Thermal Indices on the system are precise to 0.1 units.

#### 5.2.3 Acoustic Power Statement

### 5.2.3.1 Acoustic Measurement Precision and Uncertainty

When estimating accuracy of displayed numerical values, many factors are considered:

- The Probe Changeability
- The system changeability
- Changeability and accuracy of measurement
- Possible operating conditions and testing numbers needed to obtain displayed result accuracy of the diagnostic system
- Whether the display accuracy depends on specific system combination, mode combination, probe component and launch mode combination, or all of above
- Algorithm accuracy of the system software used to calculate the MI/TI
- Approximation engineering method used in real time computation

#### 5.2.3.2 Differences between Actual and Displayed MI and TI

For many assumptions used in the process of measurement and calculation, they actually are conservative. For most organizations path, high estimate is made in the measurement and calculation process of tissue exposure intensity. For example, using attenuation coefficient 0.3 dBcm<sup>-1</sup> MHz<sup>-1</sup> much lower than the actual animal tissue



attenuation coefficient, choosing conservative values of tissue characteristic.

Therefore, displayed MI and TI values should be relative information for reference, they serve to indicate to the operator whether a particular setting of the system increases or decreases the possibility of Thermal or Mechanical effect, used to help the operator be careful to use ultrasonic diagnostic system and follow the ALARA principle, these values cannot be equal to actual values.

### 5.2.3.3 Uncertainty of Measurement

Sound pressure is the most basic data of sound field measurement, and other sound field parameters can be deduced from sound pressure, so when analyzing measurement uncertainty, only take sound pressure for analysis and uncertainty of other parameters can be deduced from the sound pressure.

Measurement uncertainty mainly include repeated measurement uncertainty and the system uncertainty, the system uncertainty is an order of magnitude higher than repeated measurement uncertainty, so the main analysis is the system uncertainty. Mainly decided by the following factors:

- 1) The sensitivity of hydrophone: According to hydrophone calibration report provided by ONDA company, the maximum allowable error of sound pressure for hydrophone is ±12%
- 2) Scope: according to Agilent DSO6502A specifications, its effect on the sound pressure is ±2%
- 3) Temperature: effect of the thermocouple on sound pressure error is  $\pm 4\%$

Above all uncertainty components are not related, synthetic standard uncertainty of sound pressure is  $\pm 13\%$ 

### **5.2.3.4 Operator Control Property**

There are several system controls that the operator can use to adjust the image quality and limit the acoustic intensity. These controls are related to the techniques that an operator might use to implement ALARA and can be divided into three categories: direct, indirect, and receiver controls



#### Direct Controls

The direct control of the acoustic output of this system is adjusting voltage size. But its maximum acoustic output shouldn't be more than displayed acoustic output limit in any modes.

#### Indirect Controls

The controls that indirectly affect output are many imaging parameters. These are operating modes, frequency, focal point number/position, image depth and pulse repetition frequency (PRF).

The operating mode determines whether the ultrasound beam is scanning or non-scanning. Thermal effect is closely connected to B/M Mode, PW and Color Mode.

Acoustic attenuation of tissue is directly related to transducer frequency.

The focal point number and position is related to active aperture of transducer and beam width.

For the pulse repetition frequency (PRF), the higher the PRF, the more acoustic output power increased over a period of time.

#### Receiver Controls

The receiver control does not affect the acoustic output, including gain, dynamic range, and image processing, etc. Therefore, in the image optimization, should adjust the receiver control to optimize images firstly, the second are through direct control and indirect control.

When acquiring images, it is recommended to use the default (or as low as possible) acoustic output location and use the gain control to compensate. The default setting is commonly 70% of maximum allowed acoustic output value, which will not cause harm to the operator, and for the probe is the most effective value.

## **5.2.4 Acoustic Power Setting**

The ultrasound system has been preset the parameters for each exam mode with different probes before shipment. When the ultrasound system is powered on, a new



Animal patient is created or the application mode is changed, the system will retrieve the default settings. You can also reset the parameters.

#### **5.2.5 ALARA**

It is required to practice ALARA when using ultrasound energy. Practicing ALARA ensures that the total energy level is controlled below a low level at which bioeffects are not generated while diagnostic information is being accumulated. The total energy is controlled by output intensity and total radiation time. The output intensity necessary for examinations differs depending on the Animal patient and clinical case.

Not all examinations can be performed with an extremely low level of acoustic energy. Controlling the acoustic level at an extremely low-level leads to low-quality images or insufficient Doppler signals, adversely affecting the reliability of the diagnosis. However, the sound power which is used greater than the actual needs does not contribute to improving the quality of diagnostic information either, it will increase the risk of biological effects.

The operator must take responsibility for the safety of Animal patients.

## 5.3 Electromagnetic Compatibilities

Electromagnetic compatibilities are the abilities of the system or equipment to operate normally in the electromagnetic environment and not to radiate any electromagnetic interruptions to any other objects which are in the same environment.

This system is designed in accordance with the current EMC requirement. And the ultrasound image will degrade instantly if the system is used in the electromagnetic field environment. If the degradation of the image is found, it is recommended to inspect the operation environment to confirm the radiation source.

## 5.3.1 Electromagnetic Emission

This product is intended for use in the electromagnetic environment specified below. Users should make sure that the product is being used in such an environment.



## **Guidance and Manufacturer's Declaration – Electromagnetic Emissions**

The Wireless Probe Series Ultrasound Scanner is intended for use in the electromagnetic environment specified below. The customer or the user of Wireless Probe Series Ultrasound Scanner should assure that it is used in such environment.

environment.			
Type	Emission	Compliance	Electromagnetic environment -
Туре	Test	Compliance	guidance
RF	CISPR 11	Group1,	The Ultrasound System uses RF energy
Emission	EN 55011	Class A	for its internal functions only. Therefore,
			its RF emissions are very low and are
			not likely to cause any interference in
			nearby electronic equipment.
RF	CISPR 11	Group1,	The Ultrasound System is suitable for
Emission	EN 55011	Class A	use in all establishments, including
Harmonic	IEC 61000-	Class A	domestic establishments and those
Emission	3-2		directly connected to the public low-
	EN 61000-		voltage power supply network that
	3-2		supplies buildings used for domestic
Flicker	IEC 61000-	Pst: 1	purposes.
Emission	3-3	Plt: 0.65	
	EN 61000-	Tmax:0.5	
	3-3	dmax: 4%	
		dc: 3.3%	

< **44** / **92** > Doc # : BN\_SMV\_EOP\_01



## 5.3.2 Electromagnetic Immunity

## **Guidance and Manufacturer's Declaration – Electromagnetic Immunity**

The Wireless Probe Series Ultrasound Scanner is intended for use in the electromagnetic environment specified below. The customer or the user of the Wireless Probe Series Ultrasound Scanner should assure that it is used in such an environment.

Type	Immunity	Compliance	Electromagnetic environment -
Туре	Test	Level	guidance
Electrostatic discharge (ESD)	IEC 61000- 4-2 EN 61000- 4-2	±8 kV/Contact  ±2, ±4, ±8, ±15 kV/Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
RF Electromag netic Immunity	IEC 61000- 4-3 EN 61000- 4-3	3V/m 80 MHz - 2.7GHz 80% AM at 1kHz, 2Hz	The field strength outside the shielded location from the fixed RF transmitter as determined by electromagnetic field survey should be less than 3V/m.
Near Field Immunity for RF Wireless Communica tion Equipment	IEC 61000- 4-3 EN 61000- 4-3	Table #9 at IEC 60601-1-2: 2014	Interference may occur in the vicinity of equipment marked with the following symbol:  (((()))
Electrostatic transient burst	IEC 61000- 4-4 EN 61000- 4-4	±2 kV, 100 kHz Repeated Frequency	The quality of the mains power must be that of a typical commercial or hospital environment.
Surge	IEC 61000- 4-5 EN 61000- 4-5	(Wire to Wire) ±0.5 kV, ±1 kV	The quality of the mains power must be that of a typical commercial or hospital environment.

< **45** / **92** > Doc # : BN\_SMV\_EOP\_01



Conducted	IEC 61000-	3V	The Ultrasound System is designed for
RF	4-6	0.15 MHz - 80	use in electromagnetic environments
	EN 61000-	MHz	where radiated RF disturbances are
	4-6		controlled. Customers or users of the
		V in ISM	ultrasonic scanner series can avoid
		bands	electromagnetic interference by
		between	maintaining a minimum distance. It can
		0.15 MHz and	be used between portable and mobile RF
		80 MHz	communication equipment (transmitter)
			and ultrasonic scanner series depending
		80% AM at 1	on the maximum power of the
		kHz, 2 Hz	communication equipment as
			recommended below.
Power Frequency Magnetic Field Immunity Voltage dips	IEC 61000- 4-8 EN 61000- 4-8 IEC 61000- 4-11	30 A/m 50 Hz and 60 Hz 0 % UT: 0.5 cycle	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  Mains power quality should be that of atypical commercial or hospital
'	EN 61000-	At 0°, 45°,	environment. If the user of the Wireless
	4-11	90°, 135°, 180°, 225°, 270° and 315°	Probe Series Ultrasound Scanner requires continued operation during power mains interruptions, it is recommended that the Wireless Probe Series Ultrasound Scanner be powered from an uninterruptible power supply or a battery.
Voltage	IEC 61000-	0 % UT;	
interruption	4-11	250/300	
S	EN 61000- 4-11	cycle	

[Notice]  $\mathcal{U}_T$  is the a.c. mains voltage prior to application of the test level.



## 5.3.3 Recommended Separation Distance

The SonoMe VET is intended for use in an electromagnetic environment in which radiated RF disturbance are controlled. The customer or the user of the Wireless Probe Series Ultrasound Scanner can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SonoMe VET as recommended below, according to the maximum output power of the communications equipment.



Note If the system has generated the interference (confirmed by powering on and off the system), the qualified service personnel should solve the problem by following the steps as below

- Change the location of the affected system.
- Place SonoMe VET further away from the affected system.
- Supply power to this system in other ways other than the way used currently
- Contact Bionet's Customer Service Department ASAP

#### Warning

Use of this equipment adjacent to or stacked with other equipment may result in improper operation and should be avoided. If such use is required, this and other equipment should be observed to ensure normal operation.

#### Warning

The use of accessories, converters and cables not specified or provided by the manufacturer may increase electromagnetic emissions or decrease electromagnetic immunity of this equipment and may cause malfunction.

#### Caution

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). For use in a residential environment (which normally requires CISPR 11 Class B), this equipment may not provide adequate protection for radio frequency communications services. Users may need to take



mitigation measures such as relocating or reorienting equipment.

#### Warning

- This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
  - (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.
- This device complies with FCC radiation exposure limits set forth for an uncontrolled environment and it also complies with Part 15 of the FCC RF Rules. This equipment must not be co-located or operating in conjunction with any other antenna or transmitter. End-users and installers must be provided with antenna installation instructions and consider removing the no-collocation statement.

< 48 / 92 > Doc #: BN SMV EOP 01



# **Appendix A. Specification**

## A-1. 5C-10L-14L-5CB-10LB < Common>

Complied Standard		
EN/IEC 60601-1:2012	Medical electrical equipment Part 1: General requirements for basic safety and essential performance, BF, Non-continuous operation	
EN/IEC 60601-2-37:2015	Medical Electrical Equipment Part 2-37: Particular Requirements for the Basic Safety and Essential Performance of Ultrasonic Medical Diagnostic and Monitoring Equipment	
EN/IEC 60601-1-2:2014	Electromagnetic Compatibility Requirement and tests, Class A	
	Circumstance	
	Ambient Temperature	5°C ~ +35°C (41°F ~ 95°F)
Operation	Relative Humidity	25% ~ 80%, Non-Condensing
Operation	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3000m
	Ambient Temperature	-20°C ~ +55°C(-4°F ~ +131°F)
Storage / Shipping	Relative Humidity	25% ~ 93%, Non-Condensing
Storage / Shipping	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3000m
Safety		
Type of protection	Internal-Power Type, 3.8 VDC / 4,200 mAh	
against electric shock	Type BF Applied Parts	
Degree of protection	Non-continuous operation Operating mode: 1:2 duty cycle (ON: 5 min / OFF: 10 min)	
against electric shock	Hand-Held Type	



Degrees of protection against harmful liquid	IPX1	
Degree of safety of application	The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.	
	Wireless Parameters	
WI-FI Type	2.4G Band	
Supported Type	IEEE 802.11n(HT20)	
Modulation	IEEE 802.11n:OFDM	
Channel Number	11 Channel for 20MHz bandwidth (2,412~2,462 MHz)	
Channel Separation	5 MHz	
Wireless Power	10 mW/MHz or below	
Common Technical Parameters		
Description	- Display: Tablet PC or Mobile Phone which supports by iOS 8.0 or above, Windows System.	
Element	128	
Measurement	Length, Area, Angle, Obstetrics	
Image Frame Rate	16f/s	
Battery Operation Time	3.5 hours	
Charging Mode USB Charging		
Supported Device Brand	Apple	
Equipment name	Specific low power wireless device (Wireless device for wireless data communication system)	
Manufacturer/Country	Bionet/S.Korea	



## A-1. 5C-10L-14L-5CB-10LB < Difference among devices>

5C		
Туре	Color Doppler	
Dimension & Weight	157mm × 67mm × 29mm, 250g	
Array Type	Electronic array R60	
Probe Type	Convex	
Display Mode	B, B/M, Color, PDI, PW	
Frequency	3.2/5MHz	
Depth	90mm~305mm	
Application	Animal	
Antenna Description	Internal Antenna, 2.1dBi(Max.), AH104N2450D1	
	10L	
Туре	Color Doppler	
Dimension & Weight	157mm × 67mm × 29mm, 200g	
Array Type	Electronic array L40	
Probe Type	Linear	
Display Mode	B, B/M, Color, PDI, PW	
Frequency	7.5/10MHz	
Depth	20mm~55mm	
Application	Animal	
Antenna Description	Internal Antenna, 1.6 dBi(Max.), AF216M245001-T	
14L		
Туре	Color Doppler	
Dimension & Weight	157mm × 67mm × 29mm, 200g	
Array Type	Electronic array L25	
Probe Type	Linear	



Display Mode	B, B/M, Color, PDI, PW
Frequency	10/14MHz
Depth	20mm~55mm
Application	Peripheral Vessel, Superficial Organ
Antenna Description	Internal Antenna, 2.1dBi(Max.), AH104N2450D1
	5CB
Туре	BW
Dimension & Weight	157mm × 67mm × 29mm, 250g
Array Type	Electronic array R60
Probe Type	Convex
Display Mode	B, B/M
Frequency	3.2/5MHz
Depth	90mm~280mm
Application	Animal
Antenna Description	Internal Antenna, 1.6 dBi(Max.), AF216M245001-T
	10LB
Туре	BW
Dimension & Weight	157mm × 67mm × 29mm, 200g
Array Type	Electronic array L40
Probe Type	Linear
Display Mode	B, B/M
Frequency	7.5/10MHz
Depth	20mm~55mm
Application	Animal
Antenna Description	Internal Antenna, 1.6 dBi(Max.), AF216M245001-T



## A-2. H5C10L

Complied Standard		
EN/IEC 60601-1:2012	Medical electrical equipment Part 1: General requirements for basic safety and essential performance, BF, Non-continuous operation	
EN/IEC 60601-2-37:2015	Medical Electrical Equipment Part 2-37: Particular Requirements for the Basic Safety and Essential Performance of Ultrasonic Medical Diagnostic and Monitoring Equipment	
EN/IEC 60601-1-2:2014	Electromagnetic Compa	atibility Requirement and tests,
	Circumstance	
	Ambient Temperature	5°C ~ +35°C (41°F ~ 95°F)
Operation	Relative Humidity	25% ~ 80%, Non-Condensing
Operation	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3,000m
	Ambient Temperature	-20°C ~ +55°C(-4°F ~ +131°F)
Storage / Shipping	Relative Humidity	25% ~ 93%, Non-Condensing
Storage / Shipping	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3,000m
	Safety	
Type of protection	Internal-Power Type, 3.8 VDC / 2,800 mAh	
against electric shock	Type BF Applied Parts	
Degree of protection	Non-continuous operation Operating mode: 1:2 duty cycle (ON: 5 min / OFF: 10 min)	
against electric shock	Hand-Held Type	
Degrees of protection against harmful liquid	IPX1	
Degree of safety of	The equipment is not suitable for use in the presence of	

Doc # : BN\_SMV\_EOP\_01



application	a flammable anesthetic mixture with air, oxygen, or nitrous oxide.		
Wireless Parameters			
WI-FI Type	2.4G Band		
Supported Type	IEEE 802.11n(HT20)		
Modulation	IEEE 802.11n:OFDM		
Channel Number	11 Channel for 20MHz bandwidth (2,412~2,462 MHz)		
Channel Separation	5 MHz		
Antenna Description	Chip Antenna, 2.1dBi(Max.), AH104N2450D1		
Wireless Power	10 mW/MHz or below		
	Common Technical Parameters		
Description	- Display: Tablet PC or Mobile Phone which supports by iOS 8.0 or above, Windows System.		
Element	192		
Dimension & Weight	159mm × 69mm × 29mm, 250g		
Array Type	Electronic array R60 / L40		
Probe Type	Convex array probe, Linear array probe		
Mode	B, B/M, Color, PDI, PW		
Frequency	3.2/5MHz Convex, 7.5/10MHz Linear		
Depth	Convex 90mm~305mm, Linear 20~80mm		
Measurement	[Common] Length, Area, Angle, Velocity, HR, S/D, Depth [For Obstetrics] GA(CRL), GA(BPD), GA(GS), GA(FL), GA(HC), GA(AC), EFW(BPD), EFW(FL)		
Application	Animal		
Image Frame Rate	18f/s		
Battery Operation Time	2.5 hours		
Charging Mode	Wireless Charging		
Supported Device Brand	Apple		



Equipment name	Specific low power wireless device (Wireless device for wireless data communication system)
Manufacturer/Country Bionet/S.Korea	

< **55** / **92** > Doc # : BN\_SMV\_EOP\_01



## A-3. H5C & H10L <Common>

Complied Standard			
EN/IEC 60601-1:2012	Medical electrical equipment Part 1: General requirements for basic safety and essential performance, BF, Non-continuous operation		
EN/IEC 60601-2-37:2015	Medical Electrical Equipment Part 2-37: Particular Requirements for the Basic Safety and Essential Performance of Ultrasonic Medical Diagnostic and Monitoring Equipment		
EN/IEC 60601-1-2:2014	Electromagnetic Compa	atibility Requirement and tests,	
	Circumstance		
	Ambient Temperature	5°C ~ +35°C(41°F ~ 95°F)	
Operation	Relative Humidity	25% ~ 80%, Non-Condensing	
Operation	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3,000m	
	Ambient Temperature	-20°C ~ +55°C(-4°F ~ +131°F)	
Storage / Shipping	Relative Humidity	25% ~ 93%, Non-Condensing	
Storage / Shipping	Atmospheric Pressure, Max. Altitude	700hPa ~ 1060hPa, 3,000m	
Safety			
Type of protection	Internal-Power Type, 3.8 VDC / 5,600 mAh		
against electric shock	Type BF Applied Parts		
Degree of protection	Non-continuous operation Operating mode: 1:2 dut	on cy cycle (ON: 5 min / OFF: 10 min)	
against electric shock	Hand-Held Type		
Degrees of protection against harmful liquid	IPX1		
Degree of safety of	The equipment is not su	uitable for use in the presence of	

Doc # : BN\_SMV\_EOP\_01



application	a flammable anesthetic mixture with air, oxygen, or nitrous oxide.	
Wireless Parameters		
WI-FI Type	2.4G Band	
Supported Type	IEEE 802.11n(HT20)	
Modulation	IEEE 802.11n:OFDM	
Channel Number	11 Channel for 20MHz bandwidth (2,412~2,462 MHz)	
Channel Separation	5 MHz	
Antenna Description	Chip Antenna, 2.1dBi(Max.), AH104N2450D1	
Wireless Power	10 mW/MHz or below	
	Common Technical Parameters	
Description	Tablet PC or Mobile Phone which supports by iOS 8.0 or above, Windows System.	
Element	192	
Mode	B, B/M, Color, PDI, PW	
Measurement	[Common] Length, Area, Angle, Velocity, HR, S/D, Depth [For Obstetrics] GA(CRL), GA(BPD), GA(GS), GA(FL), GA(HC), GA(AC), EFW(BPD), EFW(FL)	
Application	Animal	
Image Frame Rate	18f/s	
Battery Operation Time	5.0 hours	
Charging Mode	USB Charging	
Supported Device Brand	Apple	
Equipment name	Specific low power wireless device (Wireless device for wireless data communication system)	
Manufacturer/Country	Bionet/S.Korea	



## A-3. H5C & H10L < Difference among devices >

	H5C
Dimension	160mm × 70mm × 23mm, 250g
Array Type	Electronic array R60
Probe Type	Convex probe
Frequency	3.2/5MHz
Depth	90mm~305mm
Application	Animal
	H10L
Dimension	160mm × 64mm × 23mm, 200g
Array Type	Electronic array L40
Probe Type	Linear probe
Frequency	7.5/10MHz
Depth	20mm~100mm
Application	Animal



# **Appendix B. Acoustic Output Data**

These data are acquired through the test report of IEC 60601-2-37.

MODE 5C B Mode

Inc	dex label	MI	7	TIS .		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.32	0	.14	0	.14	N/A
Index compo	onent value		0.14	0.14	N/A	0.14	
Acoustic Parameter	р <sub>г. "</sub> at z <sub>мі</sub> (MPa)	2.04					
S	P (mW)	STOCK AS THE	29	9.86	29	9.86	N/A
	$P_{1\times 1}$ (mW)		11	.96	1	1.96	1000000
	z <sub>s</sub> (cm)	Harling No.		N/A			
	z <sub>b</sub> (cm)				the dimension	N/A	
	ZMI (cm)	3.96	THE PARTY OF THE P				
	Z <sub>PII.</sub> 。 (cm)	3.96					
	f <sub>awf</sub> (MHz)	2.39	2	.39	2	.39	N/A
Other	prr (Hz)	1598.50	Republicani	THE REAL PROPERTY.	Library Cont	BUS EASTERN STATE	10000
Information	srr (Hz)	9.09					
	n <sub>pps</sub>	2				TO BE SHOULD	
	/pa.a at z <sub>PII.a</sub> (W/cm²)	143.14					
	/ <sub>spta.</sub> , at z <sub>PII.</sub> , or z <sub>SII.</sub> , (mW/cm <sup>2</sup> )	3.55					
	I <sub>spta</sub> at ZPII OF ZSII (mW/cm²)	7.10					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
							11.68
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< **59** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 5C B + B/M Mode

Inc	dex label	MI	7	ris		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.32	0	.21	1	.23	N/A
Index compo	onent value		B:0.14 M:N/A	B:0.14 M:0.21	N/A	B:0.14 M:1.23	
Acoustic Parameter	p <sub>r.</sub> "at z <sub>M</sub> (MPa)	2.04					
S	P (mW)		B:29.86	M:29.86	B:29.86	M:29.86	N/A
	$P_{1x1}$ (mW)	Travelle	B:1	1.96	B:	11.96	
	z <sub>s</sub> (cm)			3.06	William Land		Marin To
	z <sub>b</sub> (cm)					3.86	
	z <sub>MI</sub> (cm)	3.96					
	Z <sub>PII.</sub> (cm)	3.96	N. Texture			Jan Barre	NIE STO
	f <sub>awf</sub> (MHz)	2.39	2	.39	2	2.39	
Other	prr (Hz)	1598.50	DEN NELLE	1988 BARR			
Information	srr (Hz)	9.09					
	npps	2	BARRIS				BIN
	I <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm²)	143.14					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	250.46					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	481.92					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< **60** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 5C B + Color / B + PDI Mode

Inc	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surfac e	Below surface	
Maximum in	dex value	1.32	0	.15	(	0.15	N/A
ndex component value			B:0.11 Color: 0.04	B:0.11 Color: 0.04	N/A	B:0.11 Color: 0.04	
Acoustic Parameter	р <sub>г. a</sub> at <i>z<sub>мі</sub></i> (MPa)	2.04					
S	P (mW)		B:23.14	Color:8.31	B:23.14	Color:8.31	N/A
	$P_{1x1}$ (mW)	Year Control	B:9.27	Color:3.33	B:9.27	Color:3.33	
	z <sub>s</sub> (cm)			N/A			
	Zb (cm)				Mi Autos	N/A	de la
	z <sub>MI</sub> (cm)	3.96				Wains ba	
	ZPII. (cm)	3.96	NEW STATE				
	f <sub>awf</sub> (MHz)	B:2.39	B:2.39	Color:2.44	B:2.39	Color:2.44	N/A
Other	prr (Hz)	2957.30			Diversity of		
Information	srr (Hz)	7.04					
	n <sub>pps</sub>	2					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm²)	143.14					
	/ <sub>spta.</sub> , at z <sub>PII.</sub> , or z <sub>SII.</sub> , (mW/cm <sup>2</sup> )	14.46					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	26.87					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	B:H5.0 Color:2.5	B:H5.0 Color:2.5	B:H5.0 Color:2.5	N/A	B:H5.0 Color:2.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	2.0	2.0	2.0	N/A	2.0	N/A

< **61** / **92** > Doc # : BN\_SMV\_EOP\_01



#### MODE 5C PW Mode

Inc	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.72	0	.46	2	2.23	N/A
Index compo	onent value		N/A	0.46	N/A	2.23	the later of
Acoustic Parameter	p <sub>r.a</sub> at z <sub>MI</sub> (MPa)	1.14					
S	P (mW)		6:	5.16	6	5.16	N/A
	P <sub>1x1</sub> (mW)		N	I/A		N/A	
	z <sub>s</sub> (cm)			3.06	DEAL POR		
	z <sub>b</sub> (cm)					3.28	
	z <sub>M</sub> (cm)	3.28					
	Z <sub>PII.a</sub> (cm)	3.28					
	f <sub>awf</sub> (MHz)	2.48	2	.48	2	2.48	N/A
Other	prr (Hz)	2500.00	N. M. N. S.				
Information	srr (Hz)	N/A					
	n <sub>pps</sub>	N/A					
	/pa. at z <sub>PII. a</sub> (W/cm <sup>2</sup> )	46.12					
	/ <sub>spta, a</sub> at z <sub>PII, a</sub> or z <sub>SII, a</sub> (mW/cm <sup>2</sup> )	348.77					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	611.90					
~	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.51					
Operating	Display	40	N/A	40	N/A	40	N/A
control conditions	focus(mm) Display	30.50	I toomout	170001	(Constant)	1072	
	depth(mm)	90	N/A	90	N/A	90	N/A
	Working frequency(MHz)	2.5	N/A	2.5	N/A	2.5	N/A
	Display focus number	1	N/A	1	N/A	1	N/A
	PRF(KHz)	2.5	N/A	2.5	N/A	2.5	N/A
	SV(mm)	1	N/A	1	N/A	1	N/A

< **62** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 10L B Mode

Inc	dex label	MI	7	7S	7	TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.51	0	.19	0	.19	N/A
Index compo	onent value	Beltharine	0.19	0.19	N/A	0.19	
Acoustic Parameter	p <sub>r.</sub> at z <sub>M</sub> (MPa)	1.08					
s	P (mW)		9	.02	9	.02	N/A
	P <sub>1x1</sub> (mW)		9	.02	9	.02	
	z <sub>s</sub> (cm) z <sub>b</sub> (cm)			N/A		N/A	
	z <sub>MI</sub> (cm)	0.64					
	Z <sub>PII.</sub> (cm)	0.64	Data la	HA PYSE			1
	f <sub>awf</sub> (MHz)	4.38	4	.38	4	.38	N/A
Other	prr (Hz)	3765.30	The state of	No de contra		2012	
Information	srr (Hz)	12.99					HIVAY:
	npps	4	OF THE				
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm <sup>2</sup> )	34.71					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	54.43					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	77.80					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.18					
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2	N/A	2	N/A

< 63 / 92 > Doc #: BN\_SMV\_EOP\_01



MODE 10L B + B/M Mode

In	dex label	MI	7	ris		ΓIB	TIC
35349			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.51	0	.38	0	.47	N/A
Index compo	onent value		B:0.19 M:0.19	B:0.19 M:N/A	N/A	B:0.19 M:0.47	
Acoustic Parameter	p <sub>r.</sub> "at z <sub>M</sub> (MPa)	1.08					
s	P (mW)		B:9.02	M:9.02	B:9.02	M:9.02	N/A
	P <sub>1x1</sub> (mW)		B:	9.02	B:	9.02	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)	RANGE IN				1.16	7 333
	z <sub>MI</sub> (cm)	0.64					
	z <sub>PII.</sub> 。 (cm)	0.64					
	f <sub>awf</sub> (MHz)	4.38	4.38 4.38		.38	N/A	
Other	prr (Hz)	3765.30	S. 1912				
Information	srr (Hz)	12.99					
	//pps	4					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm <sup>2</sup> )	34.71					
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	186.55					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	238.17					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.18					
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2	N/A	2	N/A

< **64** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 10L B + Color / B + PDI Mode

Inc	dex label	MI	7	7S		TIB	TIC
			At surface	Below surface	At surfac e	Below surface	
Maximum in	dex value	0.94	0	.21		0.21	N/A
Index component value			B:0.12 Color: 0.09	B:0.12 Color: 0.09	N/A	D:0.12 Color: 0.09	
Acoustic Parameter	p <sub>r. a</sub> at z <sub>M</sub> (MPa)	2.43					
S	P (mW)		B:5.79 C	Color:2.94	B:5.79 Color:2.94		N/A
	$P_{1x1}$ (mW)		B:5.79 C	Color:2.94	B:5.79	Color:2.94	<b>1988</b>
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	0.50					S/SIL
	z <sub>PII.</sub> (cm)	0.50					
	f <sub>awf</sub> (MHz)	Color: 6.73	B:4.38 C	Color:6.73 B:4.38 Color:6.73		N/A	
Other	prr (Hz)	4000.00					
Information	srr (Hz)	8.33	MET PENY	New Ministra			
	npps	13					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm <sup>2</sup> )	215.64					
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	74.14					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm²)	115.32					
	pr at zell (MPa)	2.74					
Operating control	Display focus(mm)	6	6	6	N/A	6	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	B:H10.0 Color:6.5	B:H10.0 Color:6.5	B:H10.0 Color:6.5	N/A	B:H10.0 Color:6.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	4.0	4.0	4.0	N/A	4.0	N/A

< 65 / 92 > Doc # : BN\_SMV\_EOP\_01



MODE 10L PW Mode

Inc	dex label	MI	7	TIS .	7	TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.78	0	.21	0	.51	N/A
Index compo	onent value	EL SERVIC	0.21	N/A	N/A	0.51	To see 15
Acoustic Parameter	p <sub>r.a</sub> at z <sub>M</sub> (MPa)	2.01					
S	P (mW)		6	.46	6	.46	N/A
	P <sub>1x1</sub> (mW)		N	I/A	l N	I/A	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)	mes and head				1.22	Set Set
	z <sub>Ml</sub> (cm)	0.52				0.00	
	Z <sub>PII.</sub> (cm)	0.52		Bary Die			
	f <sub>awf</sub> (MHz)	6.61	6	.61	6	.61	N/A
Other	prr (Hz)	4000.00	EOL DEVI		MESSES CO.		
Information	srr (Hz)	N/A					TO SOLE
	n <sub>pps</sub>	N/A					1200
	/pa. at z <sub>PII.</sub> (W/cm <sup>2</sup> )	120.98					
	/ <sub>spta, a</sub> at z <sub>PII, a</sub> or z <sub>SII, a</sub> (mW/cm <sup>2</sup> )	581.58					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	730.66					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.25					
Operating control	Display focus(mm)	6	6	N/A	N/A	6	N/A
conditions	Display depth(mm)	20	20	N/A	N/A	20	N/A
	Working frequency(MHz)	6.5	6.5	N/A	N/A	6.5	N/A
	Display focus number	1	1	N/A	N/A	1	N/A
	PRF(KHz)	4.0	4.0	N/A	N/A	4.0	N/A
	SV(mm) ndicates that there	1	1	N/A	N/A	1	N/A

< **66** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 14L B Mode

In	dex label	MI	7	ris		ΓIB	TIC
			At surfac e	Below surfac e	At surface	Below surface	
Maximum in	dex value	0.41	0	.15	0	.15	N/A
Index compo	onent value		0.15	0.15	N/A	0.15	
Acoustic Parameter	ρ <sub>r.α</sub> at z <sub>M</sub> (MPa)	1.02					
S	P (mW)		5.	.17	5	.17	N/A
	$P_{1x1}$ (mW)		5.	.17	5	.17	257.03
	z <sub>s</sub> (cm)			N/A		Selection in	
	Z <sub>b</sub> (cm)				<b>SALEAR</b>	N/A	
	ZMI (cm)	0.78				REL DIS	
	Z <sub>PII.a</sub> (cm)	0.78		A STATE OF THE STA			
	f <sub>awf</sub> (MHz)	6.15	6.	15	6.15		N/A
Other	prr (Hz)	3955.90	Les de la	TE WILL BY		Feedback	
Information	srr (Hz)	12.99					
	n <sub>pps</sub>	4					
	I <sub>pa.α</sub> at z <sub>PII.α</sub> (W/cm <sup>2</sup> )	29.96					
	/ <sub>spta.α</sub> at z <sub>PII.α</sub> or z <sub>SII.α</sub> (mW/cm <sup>2</sup> )	15.73					
	I <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	22.66					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.20					
				SOME			
Operating control	Display focus(mm)	8, 15	8, 15	8, 15	N/A	8, 15	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H14.0	H14.0	H14.0	N/A	H14.0	N/A
	Display focus number	2	2	2	N/A	2	N/A



MODE 14L B + B/M Mode

Inc	dex label	MI	T	'IS		ΓIB	TIC
			At surfac e	Below surfac e	At surface	Below surface	
Maximum in	dex value	0.41	0.	15	0	.21	N/A
Index compo	ndex component value		B:0.15 M:0.14	B:0.15 M:N/A	N/A	B:0.15 M:0.21	
Acoustic Parameter	р <sub>г.ɑ</sub> at z <sub>мі</sub> (MPa)	1.02					
s	P (mW)		B:5.17	M:4.83	B:5.17	7 M:4.83	N/A
	P <sub>1x1</sub> (mW)		B:5.17	M:N/A	B:5.1	7 M:N/A	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)			Variable of		0.86	
	z <sub>M</sub> i (cm)	0.78		DETAILS.			
	ZPII.α (cm)	0.78					
	f <sub>awf</sub> (MHz)	6.15	6.	15	6	.15	N/A
Other	prr (Hz)	3955.90	I Hamile				
Information	srr (Hz)	12.99					
	n <sub>pps</sub>	B:4	NAME OF STREET	TEAST MELE			
	/ <sub>pa.α</sub> at z <sub>PII.α</sub> (W/cm²)	29.96					
	I <sub>spta.α</sub> at z <sub>PII.α</sub> or z <sub>SII.α</sub> (mW/cm <sup>2</sup> )	113.45					
	I <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	158.74					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.20					
Operating control	Display focus(mm)	8, 15	8, 15	8, 15	N/A	8, 15	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H14.0	H14.0	H14.0	N/A	H14.0	N/A
	Display focus number	2	2	2	N/A	2	N/A

< **68** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 14L B + Color / B + PDI Mode

In	dex label	MI	7	7S		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.10		27		.27	N/A
Index compo	onent value		B:0.09 Color: 0.18	B:0.09 Color: 0.18	N/A	B:0.09 Color: 0.18	
Acoustic Parameter	p <sub>r.α</sub> at z <sub>M</sub> (MPa)	2.89					
S	P (mW)	100 May 100	B:3.07 C	olor: 5.45	B:3.07 (	Color: 5.45	N/A
	P <sub>1x1</sub> (mW)		B:3.07 C	olor:5.45	B:3.07 (	Color:5.45	FIRE
	z <sub>s</sub> (cm)			N/A	Table 15	March 1985	
	z <sub>b</sub> (cm)	No less marks		I THE WASTE	With Block	N/A	
	z <sub>M</sub> (cm)	0.82	1000000	MARKET			
	Z <sub>PII.α</sub> (cm)	0.82		SEASON HI	Market Sta	THE SHAPE	1.57
	f <sub>awf</sub> (MHz)	Color:6.94	B:6.15 C	olor:6.94	B:6.15 C	Color:6.94	N/A
Other Information	prr (Hz)	Color: 8000.00					
	srr (Hz)	7.69	THE REAL PROPERTY.	hard the bar		PASSAGE OF	ing.
	n <sub>pps</sub>	B:4 Color:12					
	/ <sub>pa.α</sub> at z <sub>PII.α</sub> (W/cm²)	342.70					
	/ <sub>spta,α</sub> at z <sub>PII,α</sub> or z <sub>SII,α</sub> (mW/cm <sup>2</sup> )	32.74					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	49.19					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	3.51					
Operating control	Display focus(mm)	12	12	12	N/A	12	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	B:H14.0 Color: Fixed	B:H14. 0 Color: Fixed	B:H14. 0 Color: Fixed	N/A	B:H14.0 Color: Fixed	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	8.0	8.0	8.0	N/A	8.0	N/A

< **69** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 14L PW Mode

Index label  Maximum index value Index component value		MI	TIS		TIB		TIC
		1.06	At surfac e	Below surfac e	At surface	Below surface	
			0.33		1.48		N/A
			0.33	N/A	N/A	1.48	
Acoustic Parameter s	p <sub>r.a</sub> at z <sub>мі</sub> (MPa)	2.83					
	P (mW)		9.89		9.89		N/A
	P <sub>1x1</sub> (mW)	Taka Na Ng	N/A		N/A		
	z <sub>s</sub> (cm)		N/A				
	Zb (cm)		ANTAL TA			0.86	
	Z <sub>MI</sub> (cm)	0.82					
	Z <sub>PII.a</sub> (cm)	0.82					
	f <sub>awf</sub> (MHz)	7.09	7	.09	7.09		N/A
Other Information	prr (Hz)	4000.00		THE STATE OF	OR CALL		
	srr (Hz)	N/A					
	$n_{pps}$	N/A				The same	
	I <sub>pa.α</sub> at z <sub>PII.α</sub> (W/cm <sup>2</sup> )	287.20					
	/ <sub>spta.α</sub> at z <sub>PII.α</sub> or z <sub>SII.α</sub> (mW/cm <sup>2</sup> )	1237.00					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	1848.00					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	3.46					
					Leeds a med		
Operating control conditions	Display focus(mm)	12	12	N/A	N/A	12	N/A
	Display depth(mm)	20	20	N/A	N/A	20	N/A
	Working frequency(MHz)	Fixed	Fixed	N/A	N/A	Fixed	N/A
	Display focus number	1	1	N/A	N/A	1	N/A
	PRF(KHz)	4.0	4.0	N/A	N/A	4.0	N/A
	SV(mm)	1	1	N/A	N/A	1	N/A

See above, because of MI>1 and TI>1, 14L need to display MI and TI.



## MODE 5CB B Mode

Index label				TIS			TIB	
			МІ	Scan	Non-scan		Non-	TIC
			IVII		A <sub>aprt</sub> ≤1 cm²	A <sub>aprt</sub> >1 cm <sup>2</sup>	scan	110
Maximum index value			0.65	0.16				N/A
	$\rho_{r,\alpha}$		1.17					
	P			13.31				N/A
	min of $[P_{\alpha}(Z_s), I_{ta.\alpha}(Z_s)]$							
	Zs							
Associated acoustic	$Z_{bp}$							
parameters	$z_b$							
parameters	z at max $I_{pi.\alpha}$		4.70					
	$d_{eq}(Z_b)$							
	f <sub>awf</sub>		3.25	3.25				N/A
	Dim of	X		1.29				N/A
	A <sub>aprt</sub>	Y		1.30				N/A
	$t_d$		0.73					
	prr		1250					
	$p_r$ at max $I_{pi}$		2.10					
Other information	$d_{eq}$ at max $I_{pi}$							
	I <sub>pi.α</sub> at max MI		42.30					
	Focal	FL <sub>x</sub>						
	Length	$FL_y$						
Operating control	Depth(mm)		160	160				N/A
conditions	Freq(MHz)		3.5	3.5				N/A



MODE 5CB + B/M Mode

					TIS		TIB	
Inde	ex label		MI		Non-	scan	Non-	TIC
IIIde	sx label		IVII	Scan	A <sub>aprt</sub> ≤1 cm²	A <sub>aprt</sub> >1 cm <sup>2</sup>	scan	"
Maximum index value			0.65	0.16		0.12	0.37	N/A
	$p_{r.a}$		1.17					
	P			13.31			13.11	N/A
	min of $[P_{\alpha}(Z)]$	$[Z_s), I_{ta.\alpha}(Z_s)]$				7.68		
	Zs					2.40		
Associated acoustic	$Z_{bp}$					2.19		
parameters	$Z_b$						4.70	
parameters	z at max I <sub>pi.</sub>	α	4.70					
	$d_{eq}(Z_b)$						0.62	
	<b>f</b> <sub>awf</sub>		3.25	3.25		3.25	3.25	N/A
	Dim of	X		1.29		1.29	1.29	N/A
	A <sub>aprt</sub>	Y		1.30		1.30	1.30	N/A
	t <sub>d</sub>		0.73					
	prr		1250					
	$p_r$ at max $I_p$	i	2.10					
Other information	d <sub>eq</sub> at max I	pi					0.62	
	$I_{pi.a}$ at max	MI	42.30					
	Focal	FL <sub>x</sub>				N/A		
	Length	$FL_{y}$				5.00		
Operating control	De	pth(mm)	160	160		160	160	N/A
conditions	Fre	eq(MHz)	3.5	3.5		3.5	3.5	N/A

< **72** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE 10LB B Mode

					TIS		TIB	
Inde	ex label		MI		Non-	scan	Non-	TIC
Thu.	CX IADCI		IVII	Scan	A <sub>aprt</sub> ≤1 cm²	A <sub>aprt</sub> >1 cm <sup>2</sup>	scan	110
Maximum index value			0.75	0.036				N/A
	$p_{r.a}$		1.90					
	P			1.17				N/A
	min of $[P_{\alpha}($	$Z_s$ ), $I_{ta.\alpha}$ ( $Z_s$ )]						
	Zs							
Associated acoustic	$Z_{bp}$							
parameters	$z_b$							
parameters	z at max $I_p$	i.a	1.35					
	$d_{eq}(Z_b)$							
	f <sub>awf</sub>		6.47	6.47				N/A
	Dim of	X		0.96				N/A
	A <sub>aprt</sub>	Υ		0.44				N/A
	$t_d$		0.36					
	prr		1250					
	$p_r$ at max $I_r$	oi	2.63					
Other information	d <sub>eq</sub> at max	$I_{pi}$						
	$I_{pi.\alpha}$ at max	MI	17.30					
	Focal	FL <sub>x</sub>						
	Length	$FL_{y}$						
Operating control	De	epth(mm)	80	80				N/A
conditions	Fr	eq(MHz)	7.5	7.5				N/A



MODE 10LB B + B/M Mode

·					TIS		TIB	
Ind	ex label		М		Non-	scan	Non-	TIC
mu	ex label		IVII	Scan	A <sub>aprt</sub> ≤1 cm²	A <sub>aprt</sub> >1 cm <sup>2</sup>	scan	110
Maximum index value			0.75	0.036	0.038		0.13	N/A
	$p_{r.\alpha}$		1.90					
	P			1.17	1.22		1.22	N/A
	min of $[P_{\alpha}($	$(Z_s), I_{ta.\alpha}(Z_s)]$						
	Z <sub>S</sub>							
Associated acoustic	$Z_{bp}$							
parameters	$z_b$						1.35	
parameters	z at max $I_p$	i.α	1.35					
	$d_{eq}(Z_b)$						0.18	
	f <sub>awf</sub>		6.47	6.47	6.47		6.47	N/A
	Dim of	X		0.96	0.96		0.96	N/A
	A <sub>aprt</sub>	Y		0.44	0.44		0.44	N/A
	t <sub>d</sub>		0.36					
	prr		1250					
	$p_r$ at max I	pi	2.63					
Other information	d <sub>eq</sub> at max	I <sub>pi</sub>					0.18	
	$I_{pi.\alpha}$ at max		17.30					
	Focal	FL <sub>x</sub>			N/A			
	Length	$FL_{v}$			1.30			
Operating control	De	epth(mm)	80	80	80		80	N/A
conditions		eq(MHz)	7.5	7.5	7.5		7.5	N/A

< **74** / **92** > Doc # : BN\_SMV\_EOP\_01



These data are acquired through the test report of IEC 60601-2-37.

IEC 60601-2-37:2007 +AMD1:2015 CSV © IEC 2015

MODE H5C10L(C) B Mode

Inc	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.32	0	.14	(	).14	N/A
Index compo	onent value		0.14	0.14	N/A	0.14	
Acoustic Parameter	р <sub>г.</sub> at z <sub>мі</sub> (MPa)	2.04					
S	P (mW)		2	9.86	2	9.86	N/A
	$P_{1x1}$ (mW)		1	1.96	1	1.96	
	Z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	3.96					
	Z <sub>PII.</sub> 。 (cm)	3.96					
	f <sub>awf</sub> (MHz)	2.39	2	.39	2	2.39	N/A
Other	prr (Hz)	1598.50					
Information	srr (Hz)	9.09					
	npps	2	PATRICE				
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm²)	143.14					
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	3.55					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	7.10					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< **75** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\mathsf{H5C10L}(C)}$  B + B/M Mode

In	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.32	0	.21	1	.23	N/A
Index compo	onent value		B:0.14 M:N/A	B:0.14 M:0.21	N/A	B:0.14 M:1.23	
Acoustic Parameter	p <sub>r.</sub> at z <sub>M</sub> (MPa)	2.04					
s	P (mW)		B:29.86	M:29.86	B:29.80	6 M:29.86	N/A
	$P_{1x1}$ (mW)		B:	11.96	B:	11.96	
	z <sub>s</sub> (cm)			3.06			EN T
	z <sub>b</sub> (cm)					3.86	
	ZMI (cm)	3.96					
	ZPII.。 (cm)	3.96			Call Co		
	f <sub>awf</sub> (MHz)	2.39	2	.39	2	.39	N/A
Other	prr (Hz)	1598.50				FINANCE	
Information	srr (Hz)	9.09					
	npps	2				Mary Hon	100
	I <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm <sup>2</sup> )	143.14					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	250.46					
	I <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	481.92					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< **76** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\mathsf{H5C10L}(C)}$  B + Color / B + PDI Mode

Inc	dex label	MI	7	TIS .		TIB	TIC
			At surface	Below surface	At surfac e	Below surface	
Maximum in	dex value	1.32	0	.15	(	0.15	N/A
Index compo	onent value		B:0.11 Color: 0.04	B:0.11 Color: 0.04	N/A	B:0.11 Color: 0.04	
Acoustic Parameter	p <sub>r.</sub> at z <sub>M</sub> (MPa)	2.04					
S	P (mW)		B:23.14	Color:8.31	B:23.14	Color:8.31	N/A
	$P_{1x1}$ (mW)	francis seed	B:9.27 (	Color:3.33	B:9.27	Color:3.33	
	z <sub>s</sub> (cm)	MALINESSY		N/A			
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	3.96					
	ZPII. (cm)	3.96					
	f <sub>awf</sub> (MHz)	B:2.39	B:2.39 (	Color:2.44	B:2.39	Color:2.44	N/A
Other	prr (Hz)	2957.30			2002/60		
Information	srr (Hz)	7.04					
	npps	2					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm²)	143.14					
	/ <sub>spta.</sub> , at z <sub>PII.</sub> , or z <sub>SII.</sub> , (mW/cm²)	14.46					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	26.87					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.83					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	B:H5.0 Color:2.5	B:H5.0 Color:2.5	B:H5.0 Color:2.5	N/A	B:H5.0 Color:2.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	2.0	2.0	2.0	N/A	2.0	N/A

< 77 / 92 > Doc # : BN\_SMV\_EOP\_01



MODE <u>H5C10L(C)</u> PW Mode

Inc	dex label	MI	7	TIS .		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.72	0	.46	2	23	N/A
Index compo	onent value		N/A	0.46	N/A	2.23	
Acoustic Parameter	p <sub>r.</sub> at z <sub>M</sub> (MPa)	1.14					
S	P (mW)		65	5.16	6:	5.16	N/A
	P <sub>1x1</sub> (mW)		N	I/A	l N	I/A	
	z <sub>s</sub> (cm)			3.06		Tell Carlos	
	z <sub>b</sub> (cm)					3.28	
	z <sub>MI</sub> (cm)	3.28					
	Z <sub>PII.</sub> 。 (cm)	3.28					
	f <sub>awf</sub> (MHz)	2.48	2	.48	2	.48	N/A
Other	prr (Hz)	2500.00					Contraction of the Contraction o
Information	srr (Hz)	N/A					
	n <sub>pps</sub>	N/A					(100 pm
	I <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm <sup>2</sup> )	46.12					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	348.77					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	611.90					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.51					
Operating control	Display focus(mm)	40	N/A	40	N/A	40	N/A
conditions	Display depth(mm)	90	N/A	90	N/A	90	N/A
	Working frequency(MHz)	2.5	N/A	2.5	N/A	2.5	N/A
	Display focus number	1	N/A	1	N/A	1	N/A
	PRF(KHz)	2.5	N/A	2.5	N/A	2.5	N/A
	SV(mm)	1	N/A	1	N/A	1	N/A

< **78** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE <u>H5C10L(L)</u> B Mode

Inc	dex label	MI	7	ris		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.51	0	.19	0	).19	N/A
Index compo	onent value		0.19	0.19	N/A	0.19	
Acoustic Parameter	p <sub>r.a</sub> at z <sub>M</sub> (MPa)	1.08					
S	P (mW)		9	.02	9	0.02	N/A
	$P_{1x1}$ (mW)		9	.02	9	0.02	THE WAY
	z <sub>s</sub> (cm) z <sub>b</sub> (cm)			N/A		N/A	
	ZMI (CM)	0.64				IN/A	
	Z <sub>PII.</sub> (cm)	0.64					
	f <sub>awf</sub> (MHz)	4.38	1	.38	1	1.38	N/A
Other	prr (Hz)	3765.30		.50		r.56	14//
Information	srr (Hz)	12.99					
	npps	4					
	/pa. at z <sub>PII. a</sub> (W/cm <sup>2</sup> )	34.71					
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	54.43					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	77.80					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.18					
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2	N/A	2	N/A



MODE  $\underline{\text{H5C10L(L)}}$  B + B/M Mode

In	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.51	0	.38	0	.47	N/A
Index compo	onent value		B:0.19 M:0.19	B:0.19 M:N/A	N/A	B:0.19 M:0.47	
Acoustic Parameter	p <sub>r.a</sub> at z <sub>M</sub> (MPa)	1.08					
S	P (mW)		B:9.02	M:9.02	B:9.02	2 M:9.02	N/A
	$P_{1x1}$ (mW)	The manual	B:	9.02	B:	9.02	
	z <sub>s</sub> (cm)			N/A			A. Balla
	z <sub>b</sub> (cm)					1.16	
	z <sub>MI</sub> (cm)	0.64					
	Z <sub>PII.</sub> 。 (cm)	0.64					
	f <sub>awf</sub> (MHz)	4.38	4	.38	4	4.38	
Other	prr (Hz)	3765.30					
Information	srr (Hz)	12.99					
	npps	4					
	/ <sub>pa.</sub> , at z <sub>PII.</sub> , (W/cm²)	34.71					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	186.55					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	238.17					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.18					
	interplace of the state of the						
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2	N/A	2	N/A

< **80 / 92 >** Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\text{H5C10L(L)}}$  B + Color / B + PDI Mode

Inc	dex label	MI	7	'IS		TIB	TIC
			At surface	Below surface	At surfac e	Below surface	
Maximum in	dex value	0.94	0.	21		0.21	N/A
Index compo	onent value		B:0.12 Color: 0.09	B:0.12 Color: 0.09	N/A	B:0.12 Color: 0.09	
Acoustic Parameter	р <sub>г. «</sub> at <i>z<sub>мі</sub></i> (MPa)	2.43					
S	P (mW)		B:5.79 C	Color:2.94	B:5.79	Color:2.94	N/A
	$P_{1x1}$ (mW)		B:5.79 C	Color:2.94	B:5.79	Color:2.94	
	z <sub>s</sub> (cm)			N/A	EQLES IS		
	z <sub>b</sub> (cm)			Mary Control		N/A	
	z <sub>MI</sub> (cm)	0.50					PER HIS
	ZPII.u (cm)	0.50					
	f <sub>awf</sub> (MHz)	Color: 6.73	B:4.38 C	Color:6.73	B:4.38	Color:6.73	N/A
Other	prr (Hz)	4000.00	I North Miles				
Information	srr (Hz)	8.33					
	npps	13					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm²)	215.64					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	74.14					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	115.32					
	p <sub>r</sub> at z <sub>PII</sub> (MPa)	2.74					
Operating control	Display focus(mm)	6	6	6	N/A	6	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	B:H10.0 Color:6.5	B:H10.0 Color:6.5	B:H10.0 Color:6.5	N/A	B:H10.0 Color:6.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	4.0	4.0	4.0	N/A	4.0	N/A

< **81** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE <u>H5C10L(L)</u> PW Mode

Inc	dex label	MI	7	TIS .		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.78	0	.21	0	.51	N/A
Index compo	onent value		0.21	N/A	N/A	0.51	
Acoustic Parameter	p <sub>r.a</sub> at z <sub>M</sub> (MPa)	2.01					
S	P (mW)		6	.46	6	.46	N/A
	P <sub>1x1</sub> (mW)		N	I/A	l N	N/A	Sept 100
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					1.22	Te Miles
	z <sub>MI</sub> (cm)	0.52					
	Z <sub>PII.</sub> (cm)	0.52					H SE
	f <sub>awf</sub> (MHz)	6.61	6	.61	6	.61	N/A
Other	prr (Hz)	4000.00	ACCESS OF THE PARTY		HUE CHECK		nu Mai
Information	srr (Hz)	N/A					
	npps	N/A					
	/pa.a at z <sub>PII.a</sub> (W/cm²)	120.98					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	581.58					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	730.66					
1/2-12-12-12-12-12-12-12-12-12-12-12-12-12	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.25					
Operating control	Display focus(mm)	6	6	N/A	N/A	6	N/A
control	Display depth(mm)	20	20	N/A	N/A	20	N/A
	Working frequency(MHz)	6.5	6.5	N/A	N/A	6.5	N/A
	Display focus number	1	1	N/A	N/A	1	N/A
	PRF(KHz)	4.0	4.0	N/A	N/A	4.0	N/A
	SV(mm)	1	1	N/A	N/A	1	N/A

< **82** / **92** > Doc # : BN\_SMV\_EOP\_01



These data are acquired through the test report of IEC 60601-2-37.

IEC 60601-2-37:2007 +AMD1:2015 CSV © IEC 2015

MODE <u>H5C</u> B Mode

Inc	dex label	MI	7	7S		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.35	0	.14	0.14		N/A
Index compo	onent value		0.14	0.14	N/A	0.14	
Acoustic Parameter	p <sub>r.a</sub> at z <sub>MI</sub> (MPa)	2.08					
S	P (mW)		31	.10	3	1.10	N/A
	$P_{1x1}$ (mW)		12	2.46	1	2.46	
	z <sub>s</sub> (cm)			N/A	Medical		2002
	z <sub>b</sub> (cm)	Description				N/A	
	z <sub>Mi</sub> (cm)	3.96		State of the	28015	10/2/07/25	
	Z <sub>PII.</sub> 。 (cm)	3.96					
	f <sub>awf</sub> (MHz)	2.39	2	.39	2	2.39	N/A
Other	prr (Hz)	1598.50		PROCESSED.			
Information	srr (Hz)	9.09					
	npps	2					
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm²)	149.10			I and a		
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	3.70					
	I <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	7.40					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.89					
					EP BOSE		
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< 83 / 92 > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\mathsf{H5C}}$  B + B/M Mode

Inc	dex label	MI		TIS		TIB	TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	1.35	0	.21	1	.28	N/A
Index compo	onent value		B:0.14 M:N/A	B:0.14 M:0.21	N/A	B:0.14 M:1.28	
Acoustic Parameter	p <sub>r. a</sub> at z <sub>M</sub> (MPa)	2.08					
s	P (mW)		B:31.10	M:31.10	B:31.10	0 M:31.10	N/A
	$P_{1x1}$ (mW)		B:	12.46	B:	12.46	THE S
	z <sub>s</sub> (cm)			3.06			
	z <sub>b</sub> (cm)					3.86	
	z <sub>MI</sub> (cm)	3.96					
	z <sub>PII.</sub> (cm)	3.96					
	f <sub>awf</sub> (MHz)	2.39	2	.39	2.39		N/A
Other Information	prr (Hz)	1598.50					Ballin
	srr (Hz)	9.09					
	npps	2					
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm <sup>2</sup> )	149.10					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	260.90					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm²)	502.00					
	p <sub>r</sub> at z <sub>PII</sub> (MPa)	2.89					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	H5.0	H5.0	H5.0	N/A	H5.0	N/A
	Display focus number	1	1	1	N/A	1	N/A

< **84** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\text{H5C}}$  B + Color / B + PDI Mode

Inc	dex label	MI	7	7S		TIC	
			At surface	Below surface	At surfac e	Below surface	
Maximum in	Maximum index value		0.	.15	(	0.15	N/A
Index compo	onent value		B:0.11 Color: 0.04	B:0.11 Color: 0.04	N/A	B:0.11 Color: 0.04	
Acoustic Parameter	р <sub>г. "</sub> at <i>z<sub>мі</sub></i> (MPa)	2.08					
S	P (mW)		B:24.10	Color:8.66	B:24.10	Color:8.66	N/A
	P <sub>1×1</sub> (mW)		B:9.66 C	Color:3.47	B:9.66	Color:3.47	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	3.96				2502 500	
	ZPII. u (cm)	3.96					
	f <sub>awf</sub> (MHz)	B:2.39	B:2.39 C	Color:2.44	B:2.39 Color:2.44		N/A
Other Information	prr (Hz)	2957.30					
	srr (Hz)	7.04					
	n <sub>pps</sub>	2		WALLEY TO BE	LOYEUS 3	ES STAIN	1148.4
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm²)	149.10					
,	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	15.06					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	27.99					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.89					
Operating control	Display focus(mm)	40	40	40	N/A	40	N/A
conditions	Display depth(mm)	90	90	90	N/A	90	N/A
	Working frequency(MHz)	B:H5.0 Color:2.5	B:H5.0 Color:2.5	B:H5.0 Color:2.5	N/A	B:H5.0 Color:2.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	2.0	2.0	2.0	N/A	2.0	N/A

< **85** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE <u>H5C</u> PW Mode

Inc	Index label  Maximum index value			TIS	TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum in			0	).48	2	.32	N/A
Index compo	onent value		N/A	0.48	N/A	2.32	
Acoustic Parameter	p <sub>r. u</sub> at z <sub>Ml</sub> (MPa)	1.16					
S	P (mW)		6	7.88	6	7.88	N/A
	P <sub>1x1</sub> (mW)		1	V/A	N	J/A	VIEW S
	z <sub>s</sub> (cm)			3.06	V 200	BA BARRA	
	z <sub>b</sub> (cm)					3.28	
	z <sub>MI</sub> (cm)	3.28		Dente Story			
	Z <sub>PII.</sub> , (cm)	3.28					
	f <sub>awf</sub> (MHz)	2.48	2	.48	2	.48	N/A
Other	prr (Hz)	2500.00			The state of the		
Information	srr (Hz)	N/A					
	n <sub>pps</sub>	N/A			S ASSESSMENT		
	/ <sub>pa.</sub> at z <sub>Pil.</sub> (W/cm <sup>2</sup> )	48.04					
	/ <sub>spta, a</sub> at z <sub>PII, a</sub> or z <sub>SII, a</sub> (mW/cm <sup>2</sup> )	363.30					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	637.40					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.54					
Operating control	Display focus(mm)	40	N/A	40	N/A	40	N/A
conditions	Display depth(mm)	90	N/A	90	N/A	90	N/A
	Working frequency(MHz)	2.5	N/A	2.5	N/A	2.5	N/A
	Display focus number	1	N/A	1	N/A	1	N/A
	PRF(KHz)	2.5	N/A	2.5	N/A	2.5	N/A
	SV(mm)	1	N/A	1	N/A	1	N/A

< **86** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE <u>H10L</u> B Mode

Index label		MI		TIS	TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	Maximum index value			.20		0.20	N/A
Index compo		BEEVER STAN	0.20	0.20	N/A	0.20	Service of the servic
Acoustic Parameter	p <sub>r. a</sub> at z <sub>Ml</sub> (MPa)	1.11					
S	P (mW)		9.60		ç	0.60	N/A
	$P_{1x1}$ (mW)	English test	9	.60	g	0.60	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	0.64					
	Z <sub>PII.</sub> 。 (cm)	0.64					
	f <sub>awf</sub> (MHz)	4.38	4	.38	4	.38	N/A
Other Information	prr (Hz)	3765.30					No.
	srr (Hz)	12.99	ATTE ATTE				
	npps	4					
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm <sup>2</sup> )	36.93					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	57.90					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	82.77					
	ρ <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.22					
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2 ntended us	N/A	2	N/A

< **87** / **92** > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{\text{H10L}}$  B + B/M Mode

In	dex label	MI	7	T/S	TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum in	dex value	0.53	0.40		0	.50	N/A
Index compo	onent value		B:0.20 M:0.20	B:0.20 M:N/A	N/A	B:0.20 M:0.50	
Acoustic Parameter	p <sub>r. ₄</sub> at z <sub>M</sub> (MPa)	1.11					
S	P (mW)		B:9.60	M:9.60	B:9.60	M:9.60	N/A
	$P_{1x1}$ (mW)		B:	9.60		9.60	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)					1.16	HES
	z <sub>MI</sub> (cm)	0.64				Personal Property	
	Z <sub>PII.</sub> , (cm)	0.64					
	f <sub>awf</sub> (MHz)	4.38	4.	.38	4.38		N/A
Other Information	prr (Hz)	3765.30					
	srr (Hz)	12.99					
	n <sub>pps</sub>	4					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm²)	36.93					
	/ <sub>spta.a</sub> at z <sub>PII.a</sub> or z <sub>SII.a</sub> (mW/cm <sup>2</sup> )	198.46					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm²)	253.37					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	1.22					
	D:						
Operating control	Display focus(mm)	4, 9	4, 9	4, 9	N/A	4, 9	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	H10.0	H10.0	H10.0	N/A	H10.0	N/A
	Display focus number	2	2	2	N/A	2	N/A

< 88 / 92 > Doc # : BN\_SMV\_EOP\_01



MODE  $\underline{H10L}$  B + Color / B + PDI Mode

Inc	dex label	MI	T	IS		TIB	TIC
			At surface	Below surface	At surfac e	Below surface	
Maximum index value		0.97	0.23		(	0.23	N/A
Index compo	onent value		B:0.13 Color: 0.10	B:0.13 Color: 0.10	N/A	B:0.13 Color: 0.10	
Acoustic Parameter	p <sub>r.</sub> at z <sub>M</sub> (MPa)	2.51					
S	P (mW)		B:6.16 Color:3.13		B:6.16	Color:3.13	N/A
	P <sub>1×1</sub> (mW)		B:6.16 C	olor:3.13	B:6.16	Color:3.13	
	z <sub>s</sub> (cm)			N/A		BY KING	
	z <sub>b</sub> (cm)					N/A	
	z <sub>MI</sub> (cm)	0.50					
	Z <sub>PII.</sub> (cm)	0.50					
	f <sub>awf</sub> (MHz)	Color: 6.73	B:4.38 Color:6.73		B:4.38 Color:6.73		N/A
Other Information	prr (Hz)	4000.00		No and the same			DIFFE ST
	srr (Hz)	8.33	Madewass	Mark Settle	The Later of the L		and the same
	n <sub>pps</sub>	13					
	/ <sub>pa.</sub> at z <sub>PII.</sub> (W/cm²)	229.40					
	I <sub>spla.o</sub> at z <sub>PII.o</sub> or z <sub>SII.o</sub> (mW/cm <sup>2</sup> )	78.87					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	122.68					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.82					
Operating control	Display focus(mm)	6	6	6	N/A	6	N/A
conditions	Display depth(mm)	20	20	20	N/A	20	N/A
	Working frequency(MHz)	B:H10.0 Color:6.5	B:H10.0 Color:6.5	B:H10.0 Color:6.5	N/A	B:H10.0 Color:6.5	N/A
	Display focus number	1	1	1	N/A	1	N/A
	PRF(KHz)	4.0	4.0	4.0	N/A	4.0	N/A

< **89 / 92 >** Doc # : BN\_SMV\_EOP\_01



MODE <u>H10L</u> PW Mode

Inc	dex label	MI		TIS	TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum index value		0.80	(	0.22	(	).54	N/A
Index compo	onent value	E E E	0.22	N/A	N/A	0.54	0513161
Acoustic Parameter	p <sub>r.a</sub> at z <sub>MI</sub> (MPa)	2.07					
S	P (mW)		(	5.87	(	5.87	N/A
	$P_{1x1}$ (mW)		1	N/A		N/A	
	z <sub>s</sub> (cm)			N/A			
	z <sub>b</sub> (cm)	Meson in the		1999 14 25 15		1.22	MALTE
	z <sub>MI</sub> (cm)	0.52					
	Z <sub>PII.a</sub> (cm)	0.52					
	f <sub>awf</sub> (MHz)	6.61	(	5.61	(	5.61	N/A
Other	prr (Hz)	4000.00					
Information	srr (Hz)	N/A					
	npps	N/A					
	/ <sub>pa.a</sub> at z <sub>PII.a</sub> (W/cm <sup>2</sup> )	128.70					
	/ <sub>spta.</sub> at z <sub>PII.</sub> or z <sub>SII.</sub> (mW/cm <sup>2</sup> )	618.70					
	/ <sub>spta</sub> at z <sub>PII</sub> or z <sub>SII</sub> (mW/cm <sup>2</sup> )	777.30					
	p <sub>r.</sub> at z <sub>PII</sub> (MPa)	2.32					
Operating control	Display focus(mm)	6	6	N/A	N/A	6	N/A
conditions	Display depth(mm)	20	20	N/A	N/A	20	N/A
	Working frequency(MHz)	6.5	6.5	N/A	N/A	6.5	N/A
	Display focus number	1	1	N/A	N/A	1	N/A
	PRF(KHz)	4.0	4.0	N/A	N/A	4.0	N/A
	SV(mm)	1	1	N/A	N/A	1	N/A

< **90** / **92** > Doc # : BN\_SMV\_EOP\_01



# **PRODUCT WARRANTY**

Brand Name	SonoMeVET
Product Name	Handheld Ultrasound Scanners
Model Name	5C, 10L, 14L, 5CB, 10LB, H5C10L, H5C, H10L
Approval Number	
Approval Date	
Serial Number	
Warranty Period	
Date of Purchase	Hospital Name : Address : Name : Phone :
<b>Customer Section</b>	
Sales Agency	
Manufacturer	

- ◆ Thank you for purchasing **SonoMe VET**
- The product is manufactured and passed through strict quality control and through inspection.

< **91** / **92** > Doc # : BN\_SMV\_EOP\_01



# **Contact**



## Headquarters & International Sales & service

Bionet Co., Ltd.:

5F, 61 Digital-ro 31 gil, Guro-gu, SEOUL 08375, REPUBLIC OF KOREA

Tel: +82-2-6292-6410 / Fax: +82-2-6499-7789 / E-mail: sales@ebionet.com

Website: <a href="www.ebionet.com">www.ebionet.com</a>

#### **U.S.A sales & service representative**

Bionet America, Inc.:

2691, Dow Ave, Suite B

Tustin, CA92780 U.S.A.

Toll Free: 1-877-924-6638 / Fax: 1-714-734-1761 / E-mail: support@bionetus.com

Website: www.bionetus.com

### European sales & service representative

**Bionet Europe GmbH** 

2Li Bessemerstr. 51,

D-12103 Berlin, Germany

Tel: +49-30-240-374-52 / E-mail: bionetEU@ebionet.com

Website: <a href="https://www.ebionet.com">www.ebionet.com</a>

- End of Document -

< **92 / 92 >** Doc # : BN SMV EOP 01