









| Use | Budget | Image Quality | Model | Picture | Convex | Linear | Mode | Frequency | Element | Animal Size |
|----------|------------|---------------|------------------------|--|--------|--------|------------------------------|---------------------------------------|---------|------------------------|
| General | Minimum | Good | 5CB |  | o | | B, B/M | 3.2~5Mhz | 128 ele | Small / Medium / Large |
| General | Minimum | Good | 10LB |  | | o | B, B/M | 7.5~10Mhz | 128 ele | Small |
| Extended | Affordable | Good | 5C |  | o | | B, B/M, Color, PDI, PW | 3.2~5Mhz | 128 ele | Small / Medium / Large |
| Extended | Affordable | Good | 10L |  | | o | B, B/M, Color, PDI, PW | 7.5~10Mhz | 128 ele | Small |
| Extended | Affordable | Good | 14L |  | | o | B, B/M, Color, PDI, PW | 10~14Mhz | 128 ele | Small |
| Extended | Affordable | Better | H5C |  | o | | B, B/M, Color, PDI, PW | 3.2~5Mhz | 192 ele | Small / Medium / Large |
| Extended | Affordable | Better | H10L |  | | o | B, B/M, Color, PDI, PW | 7.5~10Mhz | 192 ele | Small |
| Extended | Reasonable | Better | H5C10L |  | o | o | B, B/M, Color, PDI, PW | Convex: 3.2~5Mhz Linear: 7.5~10Mhz | 192 ele | Small / Medium / Large |